PRODUCT CATALOG



the best version of the smile

MORELLI

morelli.com.br


45
YEARS
Orthodontics
HAS OUR NAME.

Governance and corporate continuity

TO A SOLID AND RELIABLE FUTURE

Morelli has demonstrated solidity and confidence over more than four decades of growth, focusing on long-term vision.

In 2020, it initiated a corporate succession process under the visionary leadership of its founder, Mr. Oraci Morelli. This process included thorough studies of the challenges and opportunities of succession, with expert advice on succession management. Corporate Governance and Corporate Restructuring processes have been implemented, including the creation of an experienced Advisory Board to guide decisions.

In 2023, the company established the Family Protocol and the Partners' Agreement, important steps to ensure a smooth and planned transition, ensuring that the founder's legacy endures for generations to come.

At Morelli, excellence and continuity are nonnegotiable pillars, and the commitment to clients and partners is very clear: **the future is prepared with care and dedication, and the company is ready to face it and shape it in favor of dentistry.**

That's Morelli:
Continued dedication, to always be **Morelli.**

Morelli Kits quick guide


Prescription	Model	Hook on canines	Hooks on canines and premolars	Upper canine angle	1 CASE (20 pieces) 22	1 CASE (20 pieces) 18	10 CASES (200 pieces) 22	10 CASES (200 pieces) 18	100 CASES (2000 pieces) 22	100 CASES (2000 pieces) 18	PAGE
ROTH	Self-ligating SLI	✓	✓	9°	10.14.900						11
	Self-ligating SLP	✓	✓	9°	10.13.900						11
	Self-ligating SLH	✓	✓	9°	10.23.900						11
	Ceramic SLI	✓	✓	9°	10.12.900						9
	Ceramic SLP	✓	✓	9°	10.16.900						9
	Ceramic SLH	✓	✓	9°	10.26.912						9
	Polysafira	✓	✓	9°	10.24.900						15
	Ceramic	✓	✓	9°	10.11.900						15
	Max	✓	✓	9°	10.15.903	10.15.902	10.15.933	10.15.932	10.15.923	10.15.922	20
	Max	✓	✓	13°	10.15.905	10.15.904	10.15.935	10.15.934	10.15.925	10.15.924	20
	Max (Bidimensional)	✓	✓	9°	10.15.930						20
	Nano	✓	✓	9°	10.15.906						23
	Light			13°	10.10.907	10.10.906					24
	Light			9°	10.10.915						23
	Light	✓		13°	10.10.961	10.10.960					24
	Light	✓		9°	10.10.963						24
	Light	✓	✓	13°	10.10.979	10.10.978	10.10.996				25
	Light	✓	✓	9°	10.10.981		10.10.997				25
	Light (Biofuncional)	✓	✓	13°	10.10.920						25
	Standard			13°	10.10.901	10.10.900					27
	Standard			9°	10.10.911						27
	Standard	✓		13°	10.10.951	10.10.950					28
	Standard	✓		9°	10.10.953						27
	Standard	✓	✓	13°	10.10.969						28
	Standard	✓	✓	9°	10.10.971						28
	Composite	✓	✓	13°	10.17.901		10.17.932				15
EDGEWISE	Ceramic	✓	✓		10.18.901						16
	Slim Max				10.65.905	10.65.904	10.65.935	10.65.934			22
	Slim Max	✓			10.65.907	10.65.906	10.65.937	10.65.936			23
	Slim	✓			10.65.901	10.65.900					30
	Slim				10 65 903	10 65 902					30
	Standard				10.30.901	10.30.900					
M.B.T.	Self-ligating SLI	✓	✓		10.27.900						12
	Self-ligating SLP	✓	✓		10.28.900						12
	Self-ligating SLH	✓	✓		10.23.935						12
	Ceramic	✓	✓		10.32.900						16
	Polysafira	✓	✓		10.33.900						16
	Max		✓		10.35.906		10.35.936				22
	Max	✓	✓		10.35.905		10.35.935				22
	Light				10.35.904						26
	Light	✓	✓		10.35.903						26
	Light (Bidimensional)	✓			10.35.924						26
	Standard				10.35.901						29
	Standard	✓			10.35.902						29
CAPELOZZA	Max Padrão I	✓	✓		10.19.904		10.19.934				21
	Max Padrão II	✓	✓		10.19.905		10.19.935				21
	Max Padrão III	✓	✓		10.19.906		10.19.936				21
RICKETTS	Actual	✓	✓			10.25.901					31
	Actual Especial	✓	✓			10.25.900					31
	Slim	✓	✓			10.25.905					31

Contents


FIND THE CATEGORY COLOR
AT THE PAGES EDGE

SELF-LIGATING BRACKETS 07
AESTHETIC BRACKETS 13
METAL BRACKETS 17
TUBES 33
ADHESIVES 41
ACCESSORIES 45
AESTHETIC SOLUTIONS 57
ABSOLUTE ANCHORAGE 61
BANDS 67
ARCHS 73
WIRES 81
ELASTICS 85
FUNCTIONAL ORTHOPEDICS 93
HEADCAPS AND TRACTIONERS 101
INSTRUMENTS 105
ORAL HYGIENE 117
COMPARTMENTS 121


GUIDE TO
ORTHODONTIC SPECIFICATIONS




TORQUE




SLOT



ANGLE



IN / OUT
(IN-SET / OFF-SET)




LENGTH




DISTAL OFFSET




BASE CURVATURE



DENTAL ARCH



Self-ligating Brackets



Ceramic Self-ligating


The practicality of the self-ligating technique combined **with the elegance of pure alumina ceramic**

SUPERIOR AESTHETICS

Injected in pure alumina ceramic, the translucency properties of the brackets guarantee superior aesthetics and better adaptation to tooth shades.

NITINOL CLIP

The nickel titanium alloy has elastic memory and is able to withstand mechanical deformation throughout the treatment. Due to the material's geometric memory, the Clip will always tend to return to its original position, guaranteeing a better torque reading. To complement the aesthetics of the piece, the Clip has a **matte finish**, and because it is not a coating, it does not peel and remains unchanged until the end treatment.



SIMPLE AND PRACTICAL MECHANISM

There is no need for special tools or pliers and it has a system that prevents the clip from locking due to plaque and tartar build-up.

MAXIMUM COMFORT

Considering it has small dimensions and rounded contours.

EXCELLENT ADAPTATION AND RETENTION

With optimized geometry, the base has mesiodistal and cervico-occlusal curvatures, which guarantee better adaptation to the tooth.

The base is bi-blasted and silanized, guaranteeing excellent retention, both mechanically and chemically.

Watch the Ceramic SLI presentation video

Morelli Ceramic SLI Prescription Roth with hooks on canines (9th) and premolars

3.3mm	3.3mm	3.3mm	2.9mm	3.6mm	3.6mm	2.9mm	3.3mm	3.3mm	3.3mm	Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	Angle
10.12.007	10.12.007	10.12.005	10.12.003	10.12.001	10.12.002	10.12.004	10.12.006	10.12.008	10.12.008	Replacement.022" (2 UN)*

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

3.0mm	3.0mm	3.3mm	2.7mm	2.7mm	2.7mm	2.7mm	3.3mm	3.0mm	3.0mm	Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.14.014	10.14.012	10.12.011	10.12.009	10.12.009	10.12.009	10.12.009	10.12.010	10.14.013	10.14.015	Replacement.022" (2 UN)*

Ceramic SLI
Self Ligating Interactive



1 CASE • KIT 20 PIECES
SLOT22 REF 10.12.900

Morelli Ceramic SLP Prescription Roth with hooks on canines (9th) and premolars

3.3mm	3.3mm	3.3mm	2.9mm	3.6mm	3.6mm	2.9mm	3.3mm	3.3mm	3.3mm	Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	Angle
10.16.007	10.16.007	10.16.005	10.16.003	10.16.001	10.16.002	10.16.004	10.16.006	10.16.008	10.16.008	Reposição .022" (2 UN)*

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

3.0mm	3.0mm	3.3mm	2.7mm	2.7mm	2.7mm	2.7mm	3.3mm	3.0mm	3.0mm	Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.13.014	10.13.012	10.16.011	10.16.009	10.16.009	10.16.009	10.16.009	10.16.010	10.13.013	10.13.015	Reposição .022" (2 UN)*

Ceramic SLP
Self Ligating Passive



1 CASE • KIT 20 PIECES
SLOT22 REF 10.16.900

Morelli Ceramic SLH Prescription Roth with hooks on canines (9th) and premolars

	SLP			SLI				SLP			
3.3mm	3.3mm	3.3mm	2.9mm	3.6mm	3.6mm	2.9mm	3.3mm	3.3mm	3.3mm	Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	Angle
10.16.007	10.16.007	10.16.005	10.12.003	10.12.001	10.12.002	10.12.004	10.16.006	10.16.008	10.16.008	Reposição .022" (2 UN)*

15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	

3.0mm	3.0mm	3.3mm	2.7mm	2.7mm	2.7mm	2.7mm	3.3mm	3.0mm	3.0mm	Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angulo
10.13.014	10.13.012	10.16.011	10.12.009	10.12.009	10.12.009	10.12.009	10.16.010	10.13.013	10.13.015	Reposição .022" (2 UN)*
	SLP			SLI				SLP		

Ceramic SLH
Self Ligating Hybrid



1 CASE • KIT 20 PIECES
SLOT22 REF 10.26.912

Metal lower parts 4 and 5 for better resistance to chewing.
Replacement with 2 units for ceramic parts and 5 units for metal parts.

Morelli Self-ligating line has been specially developed to offer excellent sliding mechanics, comfort for the patient and practicality for the orthodontist.

Watch the clip test video



ROUNDED DESIGN

OPTIMIZED SLOT

Slot with curvature relieve and rounded edges, reducing the effects of binding and notching. These characteristics guarantee excellent sliding mechanics.

NITINOL CLIP

The nickel-titanium alloy has elastic memory, enabling it to resist mechanical deformation throughout the entire treatment.

CATEGORY LEADER

ANTI-LOCK SYSTEM

Specially designed to prevent the Clip from jamming due to plaque and tartar buildup.

SUPERIOR RESISTANCE

Injected using the latest MIM (Metal Injection Molding) technology, its unique body is highly precise and resistant.

EXCELENTE ADAPTAÇÃO E RETENÇÃO MECÂNICA

With its optimized geometry and surface treatment, it is the most efficient base in the metal bracket market:

- Mesiodistal and cervico-occlusal curvatures, ensuring a **better adaptation;**
- Micro pins that **maximize the adhesion** area and the adhesive anchoring;
- Retention bands that **facilitate the removal of marginal excess**;
- Sandblasting for **increasing roughness.**

Morelli SLI Prescription Roth (Self-Ligating Interactive) with hooks on canines (9th) and premolars

3.0mm	3.0mm	3.0mm	2.9mm	2.9mm	2.9mm	2.9mm	3.0mm	3.0mm	3.0mm	Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	Angle
10.14.007	10.14.007	10.14.005	10.14.003	10.14.001	10.14.002	10.14.004	10.14.006	10.14.008	10.14.008	Replacement.022" (05 UN)

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

3.0mm	3.0mm	3.0mm	2.7mm	2.7mm	2.7mm	2.7mm	3.0mm	3.0mm	3.0mm	Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.14.014	10.14.012	10.14.01	10.14.009	10.14.009	10.14.009	10.14.009	10.14.011	10.14.013	10.14.015	Replacement.022" (05 UN)

SLI
Self Ligating Interactive



1 CASE • KIT 20 PIECES
SLOT22 REF 10.14.900

Morelli SLP Prescrição Roth (Autoligado passivo) com ganchos nos caninos (9°) e pré molares

3.0mm	3.0mm	3.0mm	2.9mm	2.9mm	2.9mm	2.9mm	3.0mm	3.0mm	3.0mm	Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	Angle
10.13.007	10.13.007	10.13.005	10.13.003	10.13.001	10.13.002	10.13.004	10.13.006	10.13.008	10.13.008	Replacement.022" (05 UN)

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

3.0mm	3.0mm	3.0mm	2.7mm	2.7mm	2.7mm	2.7mm	3.0mm	3.0mm	3.0mm	Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.13.014	10.13.012	10.13.010	10.13.009	10.13.009	10.13.009	10.13.009	10.13.011	10.13.013	10.13.015	Replacement.022" (05 UN)

SLP
Self Ligating Passive



1 CASE • KIT 20 PIECES
SLOT22 REF 10.13.900

Morelli SLH Prescrição Roth (Autoligado Híbrido) com ganchos nos caninos (9°) e pré molares

	SLP			SLI				SLP			
3.0mm	3.0mm	3.0mm	2.9mm	2.9mm	2.9mm	2.9mm	3.0mm	3.0mm	3.0mm	Length (Length)
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque (Torque)
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	Angle (Angle)
10.13.007	10.13.007	10.13.005	10.14.003	10.14.001	10.14.002	10.14.004	10.13.006	10.13.008	10.13.008	Replacement.022" (05 UN)

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

3.0mm	3.0mm	3.0mm	2.7mm	2.7mm	2.7mm	2.7mm	3.0mm	3.0mm	3.0mm	Length (Length)
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque (Torque)
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle (Angle)
10.13.014	10.13.012	10.13.010	10.14.009	10.14.009	10.14.009	10.14.009	10.13.011	10.13.013	10.13.015	Replacement.022" (05 UN)
	SLP			SLI				SLP		

SLH
Self Ligating Hybrid



1 CASE • KIT 20 PIECES
SLOT22 REF 10.23.900


Prescription M.B.T.

Morelli SLI Prescription M.B.T. (Self-ligating Interactive) with hooks on the canines and premolars

SLI
Self Ligating Interactive

Length	3.0mm	3.0mm	3.0mm	2.9mm	2.9mm	2.9mm	2.9mm	3.0mm	3.0mm	3.0mm
Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement.022" (5 UN)	10.27.007	10.27.007	10.27.005	10.27.003	10.27.001	10.27.002	10.27.004	10.27.006	10.27.008	10.27.008





Length	3.0mm	3.0mm	3.0mm	2.7mm	2.7mm	2.7mm	2.7mm	3.0mm	3.0mm	3.0mm
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement.022" (5 UN)	10.27.014	10.27.012	10.27.010	10.27.009	10.27.009	10.27.009	10.27.009	10.27.011	10.27.013	10.27.015


1 CASE • KIT 20 PIECES
SLOT22 REF 10.27.900

Morelli SLP Prescription M.B.T. (Self-ligating Passive) with hooks on canines and premolars

SLP
Self Ligating Passive

Length	3.0mm	3.0mm	3.0mm	2.9mm	2.9mm	2.9mm	2.9mm	3.0mm	3.0mm	3.0mm
Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement.022" (5 UN)	10.28.007	10.28.007	10.28.005	10.28.003	10.28.001	10.28.002	10.28.004	10.28.006	10.28.008	10.28.008





Length	3.0mm	3.0mm	3.0mm	2.7mm	2.7mm	2.7mm	2.7mm	3.0mm	3.0mm	3.0mm
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement.022" (5 UN)	10.28.014	10.28.012	10.28.010	10.28.009	10.28.009	10.28.009	10.28.009	10.28.011	10.28.013	10.28.015


1 CASE • KIT 20 PIECES
SLOT22 REF 10.28.900

Morelli SLH Prescription M.B.T. (Self-ligating Hybrid) with hooks on the canines and premolars

SLH
Self Ligating Hybrid

	SLP			SLI				SLP		
Length	3.0mm	3.0mm	3.0mm	2.9mm	2.9mm	2.9mm	2.9mm	3.0mm	3.0mm	3.0mm
Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement.022" (5 UN)	10.28.007	10.28.007	10.28.005	10.27.003	10.27.001	10.27.002	10.27.004	10.28.006	10.28.008	10.28.008





Length	3.0mm	3.0mm	3.0mm	2.7mm	2.7mm	2.7mm	2.7mm	3.0mm	3.0mm	3.0mm
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	2°
Replacement.022" (5 UN)	10.28.014	10.28.012	10.28.010	10.27.009	10.27.009	10.27.009	10.27.009	10.28.011	10.28.013	10.28.015
	SLP			SLI				SLP		

1 CASE • KIT 20 PIECES
SLOT22 REF 10.23.935

Aesthetic Brackets





ROTH
M.B.T.

PolySafira

- Round center hook
- Maximum polishing
- Excellent mechanical retention
- Polycrystalline sapphire
- Rounded contours
- Pure aesthetics
- Bilateral positions

More than pure aesthetics

Ceramic


- Polycrystalline, translucent and highly aesthetic
- Comfort of the rounded and polished profile
- High sliding at the edges of the slot
- Centralized ball-shaped hooks
- More comfort for the patient and less replacement code
- Application of the same bracket in bilateral positions



ROTH
M.B.T.
EDGEWISE

Morelli Polysafira Prescription Roth with hooks on canines (9th) and premolars

3.2mm	3.2mm	3.2mm	2.9mm	3.6mm	3.6mm	2.9mm	3.2mm	3.2mm	3.2mm	Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	Angle
10.24.007	10.24.007	10.24.005	10.24.003	10.24.001	10.24.002	10.24.004	10.24.006	10.24.007	10.24.007	Replacement.022" (2 UN)*



3.2mm	3.2mm	3.2mm	2.9mm	2.9mm	2.9mm	2.9mm	3.2mm	3.2mm	3.2mm	Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.24.012	10.24.011	10.24.010	10.24.008	10.24.008	10.24.008	10.24.008	10.24.009	10.24.011	10.24.012	Replacement.022" (2 UN)*




1 CASE • KIT 20 PIECES
SLOT22 REF 10.24.900

Removal tips for ceramic brackets

Morelli Ceramic Prescription Roth with hooks on canines (9th) and premolars

3.2mm	3.2mm	3.2mm	2.9mm	3.6mm	3.6mm	2.9mm	3.2mm	3.2mm	3.2mm	Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	Angle
10.11.007	10.11.007	10.11.005	10.11.003	10.11.001	10.11.002	10.11.004	10.11.006	10.11.007	10.11.007	Replacement.022" (5 UN)*



3.2mm	3.2mm	3.2mm	2.9mm	2.9mm	2.9mm	2.9mm	3.2mm	3.2mm	3.2mm	Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.11.012	10.11.011	10.11.009	10.11.008	10.11.008	10.11.008	10.11.008	10.11.010	10.11.011	10.11.012	Replacement.022" (5 UN)*




1 CASE • KIT 20 PIECES
SLOT22 REF 10.11.900

Resistance test of our ceramic pieces

Morelli Composite Prescription Roth with hooks on canines (13th) and premolars

3.2mm	3.2mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	3.2mm	3.2mm	Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°	Angle
10.17.007	10.17.007	10.17.005	10.17.003	10.17.001	10.17.002	10.17.004	10.17.006	10.17.008	10.17.008	Replacement.022" (10 UN)*



3.2mm	3.2mm	3.2mm	2.9mm	2.9mm	2.9mm	2.9mm	3.2mm	3.2mm	3.2mm	Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.17.014	10.17.012	10.17.010	10.17.009	10.17.009	10.17.009	10.17.009	10.17.011	10.17.013	10.17.015	Replacement.022" (10 UN)*



1 CASE • KIT 20 PIECES
SLOT22 REF 10.17.901

10 CASES • KIT 200 PIECES
SLOT22 REF 10.17.932

TIP: See the Sequence of Arches for Aesthetic Brackets on page. 71

Prescription ◉ M.B.T.


Morelli PolySafira Prescription M.B.T. with hooks on canines and premolars



1 CASE • KIT 20 PIECES

SLOT22 REF 10.33.900

Length	3.2mm	3.2mm	3.2mm	2.9mm	3.6mm	3.6mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement.022" (2 UN)	10.33.007	10.33.007	10.33.005	10.33.003	10.33.001	10.33.002	10.33.004	10.33.006	10.33.007	10.33.007



Length	3.2mm	3.2mm	3.2mm	2.9mm	2.9mm	2.9mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement.022" (2 UN)	10.33.013	10.33.011	10.33.009	10.33.008	10.33.008	10.33.008	10.33.008	10.33.010	10.33.012	10.33.014


Morelli Ceramic Prescription M.B.T. with hooks on canines and premolars



1 CASE • KIT 20 PIECES

SLOT22 REF 10.32.900

Length	3.2mm	3.2mm	3.2mm	2.9mm	3.6mm	3.6mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement.022" (5 UN)	10.32.007	10.32.007	10.32.005	10.32.003	10.32.001	10.32.002	10.32.004	10.32.006	10.32.007	10.32.007



Length	3.2mm	3.2mm	3.2mm	2.9mm	2.9mm	2.9mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement.022" (5 UN)	10.32.013	10.32.011	10.32.009	10.32.008	10.32.008	10.32.008	10.32.008	10.32.010	10.32.012	10.32.014

Prescrição / Prescription ◉ EDGEWISE


Morelli Ceramic Prescription Edgewise with hooks on canines and premolars



1 CASE • KIT 20 PIECES

SLOT22 REF 10.18.901

Length	3.2mm	3.2mm	3.2mm	2.9mm	3.6mm	3.6mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement.022" (10 UN)	10.18.013	10.18.013	10.18.013	10.18.012	10.18.011	10.18.011	10.18.012	10.18.013	10.18.013	10.18.013



Length	3.2mm	3.2mm	3.2mm	2.9mm	2.9mm	2.9mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement.022" (10 UN)	10.18.013	10.18.013	10.18.013	10.18.012	10.18.012	10.18.012	10.18.012	10.18.013	10.18.013	10.18.013

Metal Brackets



Brackets Line

max

Brackets from Max Line are built in a single body using the MIM (Metal Injection Molding) system, without welds, and have the strength of stainless steel.

- ✔ Patient comfort
- ✔ High-performance bonding



MORE RETENTION

Sandblasted base (35% increase in retention) and **retentive micro-pins** (16% more efficient in bonding).

CERVICO-OCCLUSAL AND MESIODISTAL CURVATURE

The bases have been improved to obtain maximum adaptation to the enamel surface and, consequently, the best anchorage.

OPENING AT THE EDGE

It allows excess resin to escape and facilitates photopolymerization.



The most complete in prescriptions

The main prescriptions are available in the Max line, all of which are produced using the MIM (Metal Injection Molding) process:

- ROTH
- CAPELOZZA
- M.B.T.
- EDGEWISE

super kit

BRACKETS
+TUBES



The best value for money
among the Morelli bracket options

2000 pieces BRACKETS max
Roth prescription
Slot 22


\+ **400 peças** TUBES FOR BONDING
All Quadrantes
Roth Prescription - Simples

- Composed of **Max Line brackets** and **Bonding Tubes**,all produced using modern MIM (Metal InjectionMolding) technology, which results in single-body parts.
- Suitable for all types of dental malocclusions, especially severe ones

Bracket Kit 100 Cases + 400 Tubes Prescription Roth – Max .022" - Can. 9th Ang. - Can / Pres. hook

2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm
-14°	-14°	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	-14°	-14°
0°	0°	0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	0°	0°
14°	14°											14°	14°
20.11.211	20.11.211	10.15.413	10.15.413	10.15.409	10.15.402	10.15.400	10.15.401	10.15.403	10.15.411	10.15.415	10.15.415	20.11.212	20.11.212

Length · Torque · Angle · Offset · Replacement.022" (10 UN)



100 CASES • 2400 PIECES
SLOT22 REF 10.15.723


2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm
-30°	-30°	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	-30°	-30°
0°	0°	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	0°	0°
4°	4°											4°	4°
20.11.221	20.11.221	10.15.426	10.15.422	10.15.418	10.15.416	10.15.416	10.15.416	10.15.416	10.15.420	10.15.424	10.15.428	20.11.222	20.11.222

Length · Torque · Angle · Offset · Replacement.022" (10 UN)

Bracket Kit 100 Cases + 400 Tubes Prescription Roth – Max .022" - Can. 13th Ang. - Can / Pres. hook

2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm
-14°	-14°	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	-14°	-14°
0°	0°	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°	0°	0°
14°	14°											14°	14°
20.11.211	20.11.211	10.15.413	10.15.413	10.15.405	10.15.402	10.15.400	10.15.401	10.15.403	10.15.407	10.15.415	10.15.415	20.11.212	20.11.212

Length · Torque · Angle · Offset · Replacement.022" (10 UN)



100 CASES • 2400 PIECES
SLOT22 REF 10.15.725

2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm
-30°	-30°	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	-30°	-30°
0°	0°	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	0°	0°
4°	4°											4°	4°
20.11.221	20.11.221	10.15.426	10.15.422	10.15.418	10.15.416	10.15.416	10.15.416	10.15.416	10.15.420	10.15.424	10.15.428	20.11.222	20.11.222

Length · Torque · Angle · Offset · Replacement.022" (10 UN)


Prescription ◉ ROTH

max

Morelli Max Prescription Roth with hooks on canines (13th) and premolars



Length	2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
Replacement .022" (10 UN)	10.15.413	10.15.413	10.15.405	10.15.402	10.15.400	10.15.401	10.15.403	10.15.407	10.15.415	10.10.415
Replacement .018" (10UN)	10.15.313	10.15.313	10.15.305	10.15.302	10.15.300	10.15.301	10.15.303	10.15.307	10.15.315	10.15.315



Length	2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement .022" (10 UN)	10.15.426	10.15.422	10.15.418	10.15.416	10.15.416	10.15.416	10.15.416	10.15.420	10.15.424	10.15.428
Replacement .018" (10 UN)	10.15.326	10.15.322	10.15.318	10.15.316	10.15.316	10.15.316	10.15.316	10.15.320	10.15.324	10.15.328


1 CASE • KIT 20 PIECES	
SLOT22	REF 10.15.905
SLOT18	REF 10.15.904
10 CASES • 200 PIECES	
SLOT22	REF 10.15.935
SLOT18	REF 10.15.934
100 CASES • 2000 PIECES	
SLOT22	REF 10.15.925
SLOT18	REF 10.15.924

max

Morelli Max Prescription Roth with hooks on canines (9th) and premolars



Length	2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°
Replacement .022" (10 UN)	10.15.413	10.15.413	10.15.409	10.15.402	10.15.400	10.15.401	10.15.403	10.15.411	10.15.415	10.10.415
Replacement .018" (10 UN)	10.15.313	10.15.313	10.15.309	10.15.302	10.15.300	10.15.301	10.15.303	10.15.311	10.15.315	10.15.315



Length	2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement .022" (10 UN)	10.15.426	10.15.422	10.15.418	10.15.416	10.15.416	10.15.416	10.15.416	10.15.420	10.15.424	10.15.428
Replacement .018" (10 UN)	10.15.326	10.15.322	10.15.318	10.15.316	10.15.316	10.15.316	10.15.316	10.15.320	10.15.324	10.15.328


1 CASe • KIT 20 PIECES	
SLOT22	REF 10.15.903
SLOT18	REF 10.15.902
10 CASES • 200 PIECES	
SLOT22	REF 10.15.933
SLOT18	REF 10.15.932
100 CASES • 2000 PIECES	
SLOT22	REF 10.15.923
SLOT18	REF 10.15.922

max

Morelli Max Two-dimensional Prescription Roth with hooks on canines (9th) and premolars



Length	2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°
Replacement .022" (10 UN)	10.15.413	10.15.413	10.15.409	10.15.302	10.15.300	10.15.301	10.15.303	10.15.411	10.15.415	10.10.415




Length	2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement" .022" (10 UN)	10.15.426	10.15.422	10.15.418	10.15.316	10.15.316	10.15.316	10.15.316	10.15.420	10.15.424	10.15.428

1 CASE • KIT 20 PIECES	
SLOT22 SLOT18	REF 10.15.930

Prescription CAPELOZZA

Morelli Max Prescription Capelozza Standard I with hooks on canines and premolars

2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm	Length
-7°	-7°	-5°	+3°	+7°	+7°	+3°	-5°	-7°	-7°	Torque
0°	0°	+8°	+9°	+5°	+5°	+9°	+8°	0°	0°	Angle
10.19.107	10.19.107	10.19.105	10.19.103	10.19.101	10.19.102	10.19.104	10.19.106	10.19.108	10.19.108	Replacement.022" (10 UN)



2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm	Length
-22°	-17°	-11°	-1°	-1°	-1°	-1°	-11°	-17°	-22°	Torque
+2°	+2°	+5°	+2°	+2°	+2°	+2°	+5°	+2°	+2°	Angle
10.19.115	10.19.113	10.19.111	10.19.109	10.19.109	10.19.110	10.19.110	10.19.112	10.19.114	10.19.116	Replacement.022" (10 UN)

max




1 CASE • KIT 20 PIECES
SLOT22 REF 10.19.904

10 CASES • 200 PIECES
SLOT22 REF 10.19.934

Morelli Max Prescription Capelozza Standard II with hooks on canines and premolars

2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm	Length
-7°	-7°	-5°	+3°	+7°	+7°	+3°	-5°	-7°	-7°	Torque
0°	0°	+5°	+9°	+5°	+5°	+9°	+5°	0°	0°	Angle
10.19.107	10.19.107	10.19.117	10.19.103	10.19.101	10.19.102	10.19.104	10.19.118	10.19.108	10.19.108	Replacement.022" (10 UN)



2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm	Length
-22°	-17°	-11°	-1°	-1°	+4°	+4°	-11°	-17°	-22°	Torque
+2°	+2°	+5°	+2°	+2°	+2°	+2°	+5°	+2°	+2°	Angle
10.19.115	10.19.113	10.19.111	10.19.119	10.19.119	10.19.120	10.19.120	10.19.112	10.19.114	10.19.116	Replacement.022" (10 UN)

max




1 CASE • KIT 20 PIECES
SLOT22 REF 10.19.905

10 CASES • 200 PIECES
SLOT22 REF 10.19.935

Morelli Max Prescription Capelozza Standard III with hooks on canines and premolars

2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm	Length
-7°	-7°	-2°	+10°	+14°	+14°	+10°	-2°	-7°	-7°	Torque
0°	0°	+11°	+9°	+5°	+5°	+9°	+11°	0°	0°	Angle
10.19.107	10.19.107	10.19.125	10.19.123	10.19.121	10.19.122	10.19.124	10.19.126	10.19.108	10.19.108	Replacement.022" (10 UN)



2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm	Length
-22°	-17°	-11°	-6°	-6°	-6°	-6°	-11°	-17°	-22°	Torque
+2°	+2°	0°	0°	0°	0°	0°	0°	+2°	+2°	Angle
10.19.115	10.19.113	10.19.128	10.19.127	10.19.127	10.19.127	10.19.127	10.19.129	10.19.114	10.19.116	Replacement.022" (10 UN)

max



1 CASE • KIT 20 PIECES
SLOT22 REF 10.19.906

10 CASES • 200 PIECES
SLOT22 REF 10.19.936

Prescription M.B.T.

Morelli Max Prescription M.B.T. with hooks on the premolars


max



1 CASe • KIT 20 PIECES
SLOT22 REF 10.35.906

10 CASES • 200 PIECES
SLOT22 REF 10.35.936

Length	2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm
Torque	-7°	-7°	-7°	+10°	+17°	+17°	+10°	-7°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement.022" (10 UN)	10.35.109	10.35.109	10.35.105	10.35.103	10.35.101	10.35.102	10.35.104	10.35.106	10.35.110	10.35.110



Length	2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm
Torque	-17°	-12°	-6°	-6°	-6°	-6°	-6°	-6°	-12°	-17°
Angle	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement.022" (10 UN)	10.35.118	10.35.116	10.35.112	10.35.111	10.35.111	10.35.111	10.35.111	10.35.113	10.35.117	10.35.119

Morelli Max Prescription M.B.T. with hooks on canines and premolars


max



1 CASE • KIT 20 PIECES
SLOT22 REF 10.35.905

10 CASES • 200 PIECES
SLOT22 REF 10.35.935

Length	2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm
Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement.022" (10 UN)	10.35.109	10.35.109	10.35.107	10.35.103	10.35.101	10.35.102	10.35.104	10.35.108	10.35.110	10.35.110



Length	2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement .022" (10 UN)	10.35.118	10.35.116	10.35.114	10.35.111	10.35.111	10.35.111	10.35.111	10.35.115	10.35.117	10.35.119

Prescription EDGEWISE

Morelli Slim Max Prescription Edgewise


Slim max



1 CASE • KIT 20 PIECES
SLOT22 REF 10.65.905
SLOT18 REF 10.65.904

10 CASES • 200 PIECES
SLOT22 REF 10.65.935
SLOT18 REF 10.65.934

Length	2.5mm	2.5mm	2.5mm	2.2mm	3.0mm	3.0mm	2.2mm	2.5mm	2.5mm	2.5mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement .022" (10 UN)	10.65.403	10.65.403	10.65.403	10.65.402	10.65.401	10.65.401	10.65.402	10.65.403	10.65.403	10.65.403
Replacement .018" (10 UN)	10.65.303	10.65.303	10.65.303	10.65.302	10.65.301	10.65.301	10.65.302	10.65.303	10.65.303	10.65.303



Length	2.5mm	2.5mm	2.5mm	2.2mm	2.2mm	2.2mm	2.2mm	2.5mm	2.5mm	2.5mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement .022" (10 UN)	10.65.403	10.65.403	10.65.403	10.65.402	10.65.402	10.65.402	10.65.402	10.65.403	10.65.403	10.65.403
Replacement .018" (10 UN)	10.65.303	10.65.303	10.65.303	10.65.302	10.65.302	10.65.302	10.65.302	10.65.303	10.65.303	10.65.303

Prescription EDGEWISE

Morelli Slim Max Prescription Edgewise with hooks on canine

2.5mm	2.5mm	2.5mm	2.2mm	3.0mm	3.0mm	2.2mm	2.5mm	2.5mm	2.5mm	
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	Torque
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	Angle
10.65.403	10.65.403	10.65.404	10.65.402	10.65.401	10.65.401	10.65.402	10.65.405	10.65.403	10.65.403	Replacement.022" (10 UN)
10.65.303	10.65.303	10.65.304	10.65.302	10.65.301	10.65.301	10.65.302	10.65.305	10.65.303	10.65.303	Replacement.018" (10 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	

2.5mm	2.5mm	2.5mm	2.2mm	2.2mm	2.2mm	2.2mm	2.5mm	2.5mm	2.5mm	
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	Torque
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	Angle
10.65.403	10.65.403	10.65.405	10.65.402	10.65.402	10.65.402	10.65.402	10.65.404	10.65.403	10.65.403	Replacement.022" (10 UN)
10.65.303	10.65.303	10.65.305	10.65.302	10.65.302	10.65.302	10.65.302	10.65.304	10.65.303	10.65.303	Replacement.018" (10 UN)

Slim max



1 CASE • KIT 20 PIECES

SLOT22	REF 10.65.907
SLOT18	REF 10.65.906

10 CASES • 200 PIECES

SLOT22	REF 10.65.937
SLOT18	REF 10.65.936

Prescription ROTH

Morelli Nano Nickel Free Prescription Roth with hooks on canines and premolars

2.5mm	2.5mm	2.6mm	2.5mm	3.2mm	3.2mm	2.5mm	2.6mm	2.5mm	2.9mm	
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	Angle
10.15.509	10.15.509	10.15.505	10.15.503	10.15.501	10.15.502	10.15.504	10.15.506	10.15.510	10.15.510	Replacement.022" (10 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	

2.5mm	2.5mm	2.6mm	2.2mm	2.2mm	2.2mm	2.2mm	2.6mm	2.5mm	2.5mm	
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.15.516	10.15.514	10.15.512	10.15.511	10.15.511	10.15.511	10.15.511	10.15.513	10.15.515	10.15.517	Replacement.022" (10 UN)

nano NICKELFREE



1 CASE • KIT 20 PIECES

SLOT22	REF 10.15.906

Morelli Light Prescription Roth (Canines 9th)

2.5mm	2.5mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm	
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	Angle
10.10.408	10.10.408	10.10.406	10.10.402	10.10.400	10.10.401	10.10.403	10.10.407	10.10.408	10.10.408	Replacement.022" (10 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	

2.7mm	2.7mm	2.7mm	2.3mm	2.3mm	2.3mm	2.3mm	2.7mm	2.7mm	2.7mm	
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.10.416	10.10.414	10.10.412	10.10.410	10.10.410	10.10.410	10.10.410	10.10.413	10.10.414	10.10.416	Replacement.022" (10 UN)

Light



1 CASE • KIT 20 PIECES

SLOT22	REF 10.10.915

Prescription ◉ ROTH


Morelli Light Prescription Roth (Canines 13th)

Light



1 CASE • KIT 20 PIECES
SLOT22 REF 10.10.907
SLOT18 REF 10.10.906

Length	2.5mm	2.5mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
Replacement.022"(10 UN)	10.10.408	10.10.408	10.10.404	10.10.402	10.10.400	10.10.401	10.10.403	10.10.405	10.10.408	10.10.408
Replacement.018"(10 UN)	10.10.308	10.10.308	10.10.304	10.10.302	10.10.300	10.10.301	10.10.303	10.10.305	10.10.308	10.10.308



Length	2.7mm	2.7mm	2.7mm	2.3mm	2.3mm	2.3mm	2.3mm	2.7mm	2.7mm	2.7mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement.022"(10 UN)	10.10.416	10.10.414	10.10.412	10.10.410	10.10.410	10.10.410	10.10.410	10.10.413	10.10.414	10.10.416
Replacement.018"(10 UN)	10.10.316	10.10.314	10.10.312	10.10.310	10.10.310	10.10.310	10.10.310	10.10.313	10.10.314	10.10.316


Morelli Light Prescription Roth with hooks on the canines (9th)

Light



1 CASE • KIT 20 PIECES
SLOT22 REF 10.10.963

Length	2.5mm	2.5mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°
Replacement.022"(10 UN)	10.10.408	10.10.408	10.10.422	10.10.402	10.10.400	10.10.401	10.10.403	10.10.423	10.10.408	10.10.408



Length	2.7mm	2.7mm	2.7mm	2.3mm	2.3mm	2.3mm	2.3mm	2.7mm	2.7mm	2.7mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement.022"(10 UN)	10.10.416	10.10.414	10.10.426	10.10.410	10.10.410	10.10.410	10.10.410	10.10.427	10.10.414	10.10.416


Morelli Light Prescription Roth with hooks on the canines (13)

Light



1 CASE • KIT 20 PIECES
SLOT22 REF 10.10.961
SLOT18 REF 10.10.960


Length	2.5mm	2.5mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
Replacement.022"(10 UN)	10.10.408	10.10.408	10.10.420	10.10.402	10.10.400	10.10.401	10.10.403	10.10.421	10.10.408	10.10.408
Replacement.018"(10 UN)	10.10.308	10.10.308	10.10.320	10.10.302	10.10.300	10.10.301	10.10.303	10.10.321	10.10.308	10.10.308



Length	2.7mm	2.7mm	2.7mm	2.3mm	2.3mm	2.3mm	2.3mm	2.7mm	2.7mm	2.7mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement.022"(10 UN)	10.10.416	10.10.414	10.10.426	10.10.410	10.10.410	10.10.410	10.10.410	10.10.427	10.10.414	10.10.416
Replacement.018"(10 UN)	10.10.316	10.10.314	10.10.326	10.10.310	10.10.310	10.10.310	10.10.310	10.10.327	10.10.314	10.10.316

Morelli Light Prescription Roth with hooks on canines (9th) and premolars

2.9mm	2.9mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.9mm	2.9mm	
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+9°	+8°	+5°	+5°	+9°	+9°	0°	0°	Angle
10.10.424	10.10.424	10.10.422	10.10.402	10.10.400	10.10.401	10.10.403	10.10.423	10.10.425	10.10.425	Replacement.022" (10 un)



2.9mm	2.9mm	2.7mm	2.3mm	2.3mm	2.3mm	2.3mm	2.7mm	2.9mm	2.9mm	
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.10.430	10.10.428	10.10.426	10.10.410	10.10.410	10.10.410	10.10.410	10.10.427	10.10.429	10.10.431	Replacement.022" (10 un)

Light




1 CASE • KIT 20 PIECES
SLOT22 REF 10.10.981

10 CASES • 200 PIECES
SLOT22 REF 10.10.997

Morelli Light Prescription Roth with hooks on canines (13th) and premolars

2.9mm	2.9mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.9mm	2.9mm	
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°	Angle
10.10.424	10.10.424	10.10.420	10.10.402	10.10.400	10.10.401	10.10.403	10.10.421	10.10.425	10.10.425	Replacement.022" (10 un)
10.10.324	10.10.324	10.10.320	10.10.302	10.10.300	10.10.301	10.10.303	10.10.321	10.10.325	10.10.325	Replacement.018" (10 un)



2.9mm	2.9mm	2.7mm	2.3mm	2.3mm	2.3mm	2.3mm	2.7mm	2.9mm	2.9mm	
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.10.430	10.10.428	10.10.426	10.10.410	10.10.410	10.10.410	10.10.410	10.10.427	10.10.429	10.10.431	Replacement.022" (10 un)
10.10.330	10.10.328	10.10.326	10.10.310	10.10.310	10.10.310	10.10.310	10.10.327	10.10.329	10.10.331	Replacement.018" (10 un)

Light




1 CASE • KIT 20 PIECES
SLOT22 REF 10.10.979
SLOT18 REF 10.10.978

10 CASES • 200 PIECES
SLOT22 REF 10.10.996

Morelli Light Roth Biofunctional Prescription Class III with hooks on canines (13th) and premolars

2.9mm	2.9mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.9mm	2.9mm	
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°	Angle
10.10.424	10.10.424	10.10.420	10.10.436	10.10.434	10.10.435	10.10.437	10.10.421	10.10.425	10.10.425	Replacement.022" (10 un)



2.9mm	2.9mm	2.7mm	2.2mm	2.2mm	2.2mm	2.2mm	2.7mm	2.9mm	2.9mm	
-22°	-17°	-11°	+10°	+10°	+10°	+10°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.10.430	10.10.428	10.10.426	10.10.432	10.10.438	10.10.438	10.10.433	10.10.427	10.10.429	10.10.431	Replacement.022" (10 un)

Light



1 CASE • KIT 20 PIECES
SLOT22 REF 10.10.920

Prescription M.B.T.

Morelli Light HP Two-dimensional Prescription M.B.T. with hooks on canine


Light HP



1 CASE • KIT 20 PIECES

SLOT22 SLOT18 REF 10.35.924

Length	2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm
Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement.022"(10 UN)	10.35.412	10.35.412	10.35.417	10.35.403	10.35.401	10.35.402	10.35.404	10.35.418	10.35.412	10.35.412



Length	2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle	0°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement.022"(10 UN)	10.35.415	10.35.413	10.35.420	10.35.419	10.35.419	10.35.419	10.35.419	10.35.421	10.35.414	10.35.416

Morelli Light Prescrição M.B.T.


Light



1 CASE • KIT 20 PIECES

SLOT22 REF 10.35.904

Length	2.9mm	2.9mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.9mm	2.9mm
Torque	-7°	-7°	-7°	+10°	+17°	+17°	+10°	-7°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement.022"(10 UN)	10.35.007	10.35.007	10.35.030	10.35.103	10.35.101	10.35.102	10.35.104	10.35.031	10.35.007	10.35.007



Length	2.9mm	3.2mm	3.2mm	2.3mm	2.3mm	2.3mm	2.3mm	3.2mm	3.2mm	2.9mm
Torque	-17°	-12°	-6°	-6°	-6°	-6°	-6°	-6°	-12°	-17°
Angle	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement.022"(10 UN)	10.35.015	10.35.013	10.35.032	10.35.111	10.35.111	10.35.111	10.35.111	10.35.033	10.35.014	10.35.016

Morelli Light Prescription M.B.T. with hooks on canines and premolars


Light



1 CASE • KIT 20 PIECES

SLOT22 REF 10.35.903

Length	2.9mm	2;9mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.9mm	2;9mm
Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement.022"(10 UN)	10.35.034	10.35.034	10.35.025	10.35.103	10.35.101	10.35.102	10.35.104	10.35.026	10.35.035	10.35.035



Length	2.9mm	2.9mm	2.7mm	2.3mm	2.3mm	2.3mm	2.3mm	2.7mm	2.9mm	2.9mm
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement.022"(10 UN)	10.35.038	10.35.036	10.35.028	10.35.111	10.35.111	10.35.111	10.35.111	10.35.029	10.35.037	10.35.039

Morelli Standard Prescription Roth (Canines 9th)

STANDARD

2.5mm	2.5mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.5mm	2.5mm	
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	Angle
10.10.208	10.10.208	10.10.206	10.10.202	10.10.200	10.10.201	10.10.203	10.10.207	10.10.208	10.10.208	Replacement.022" (10 UN)

15 14 13 12 11 21 22 23 24 25



45 44 43 42 41 31 32 33 34 35

2.7mm	2.7mm	3.2mm	2.3mm	2.3mm	2.3mm	2.3mm	3.2mm	2.7mm	2.7mm	
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.10.216	10.10.214	10.10.212	10.10.410	10.10.410	10.10.410	10.10.410	10.10.213	10.10.214	10.10.216	Replacement.022" (10 UN)



1 CASE • KIT 20 PIECES

SLOT22 REF 10.10.911

Morelli Standard Prescription Roth (Canines 13th)

STANDARD

2.5mm	2.5mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.5mm	2.5mm	
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°	Angle
10.10.208	10.10.208	10.10.204	10.10.202	10.10.200	10.10.201	10.10.203	10.10.205	10.10.208	10.10.208	Replacement.022" (10 UN)
10.10.108	10.10.108	10.10.104	10.10.102	10.10.100	10.10.101	10.10.103	10.10.105	10.10.108	10.10.108	Replacement.018" (10 UN)

15 14 13 12 11 21 22 23 24 25



45 44 43 42 41 31 32 33 34 35

2.7mm	2.7mm	3.2mm	2.3mm	2.3mm	2.3mm	2.3mm	3.2mm	2.7mm	2.7mm	
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.10.216	10.10.214	10.10.212	10.10.410	10.10.410	10.10.410	10.10.410	10.10.213	10.10.214	10.10.216	Replacement.022" (10 UN)
10.10.116	10.10.114	10.10.112	10.10.310	10.10.310	10.10.310	10.10.310	10.10.113	10.10.114	10.10.116	Replacement.018" (10 UN)



1 CASE • KIT 20 PIECES

SLOT22 REF 10.10.901

SLOT18 REF 10.10.900

Morelli Standard Prescription Roth with hooks on the canines (9th)

STANDARD

2.5mm	2.5mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.5mm	2.5mm	
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	Angle
10.10.208	10.10.208	10.10.222	10.10.202	10.10.200	10.10.201	10.10.203	10.10.223	10.10.208	10.10.208	Replacement.022" (10 UN)

15 14 13 12 11 21 22 23 24 25



45 44 43 42 41 31 32 33 34 35

2.7mm	2.7mm	3.2mm	2.3mm	2.3mm	2.3mm	2.3mm	3.2mm	2.7mm	2.7mm	
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Angle
10.10.216	10.10.214	10.10.226	10.10.410	10.10.410	10.10.410	10.10.410	10.10.227	10.10.214	10.10.216	Replacement.022" (10 UN)



1 CASE • KIT 20 PIECES

SLOT22 REF 10.10.953

Prescription ◉ ROTH

Morelli Standard Prescription Roth with hooks on the canines (13th)

STANDARD




1 CASE • KIT 20 PIECES

SLOT22 REF 10.10.951

SLOT18 REF 10.10.950

Length	2.5mm	2.5mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.5mm	2.5mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
Replacement.022"(10 UN)	10.10.208	10.10.208	10.10.220	10.10.202	10.10.200	10.10.201	10.10.203	10.10.221	10.10.208	10.10.208
Replacement.018"(10 UN)	10.10.108	10.10.108	10.10.120	10.10.102	10.10.100	10.10.101	10.10.103	10.10.121	10.10.108	10.10.108



Length	2.7mm	2.7mm	3.2mm	2.3mm	2.3mm	2.3mm	2.3mm	3.2mm	2.7mm	2.7mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement.022"(10 UN)	10.10.216	10.10.214	10.10.226	10.10.410	10.10.410	10.10.410	10.10.410	10.10.227	10.10.214	10.10.216
Replacement.018"(10 UN)	10.10.116	10.10.114	10.10.126	10.10.310	10.10.310	10.10.310	10.10.310	10.10.127	10.10.114	10.10.116

Morelli Standard Prescription Roth with hooks on canines (9th) and premolars


STANDARD



1 CASE • KIT 20 PIECES

SLOT22 REF 10.10.971

Length	2.9mm	2.9mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.9mm	2.9mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°
Replacement.022"(10 UN)	10.10.224	10.10.224	10.10.222	10.10.202	10.10.200	10.10.201	10.10.203	10.10.223	10.10.225	10.10.225



Length	2.9mm	2.9mm	3.2mm	2.3mm	2.3mm	2.3mm	2.3mm	3.2mm	2.9mm	2.9mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement.022"(10 UN)	10.10.230	10.10.228	10.10.226	10.10.410	10.10.410	10.10.410	10.10.410	10.10.227	10.10.229	10.10.231

Morelli Standard Prescription Roth with hooks on canines (13th) and premolars


STANDARD



1 CASE • KIT 20 PIECES

SLOT22 REF 10.10.969

Length	2.9mm	2.9mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.9mm	2.9mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
Replacement.022"(10 UN)	10.10.224	10.10.224	10.10.220	10.10.202	10.10.200	10.10.201	10.10.203	10.10.221	10.10.225	10.10.225




Length	2,9mm	2.9mm	3.2mm	2.3mm	2.3mm	2.3mm	2.3mm	3.2mm	2.9mm	2.9mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	+2°	0°
Replacement.022"(10 UN)	10.10.230	10.10.228	10.10.226	10.10.410	10.10.410	10.10.410	10.10.410	10.10.227	10.10.229	10.10.231

Morelli Standard Prescription M.B.T.

STANDARD

2.9mm	2.9mm	2.9mm	2.9mm	3.9mm	3.9mm	2.9mm	2.9mm	2.9mm	2.9mm	Length
-7°	-7°	-7°	+10°	+17°	+17°	+10°	-7°	-7°	-7°	Torque
0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°	Angle
10.35.007	10.35.007	10.35.005	10.35.003	10.35.001	10.35.002	10.35.004	10.35.006	10.35.007	10.35.007	Replacement.022" (10 UN)



2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm	Length
-17°	-12°	-6°	-6°	-6°	-6°	-6°	-6°	-12°	-17°	Torque
+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°	Angle
10.35.015	10.35.013	10.35.011	10.35.111	10.35.111	10.35.111	10.35.111	10.35.012	10.35.014	10.35.016	Replacement.022" (10 UN)




1 CASE • KIT 20 PIECES

SLOT22 REF 10.35.901

Morelli Standard Prescription M.B.T. with hooks on the canines

STANDARD

2.9mm	2.9mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.9mm	2.9mm	Length
-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°	Torque
0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°	Angle
10.35.007	10.35.007	10.35.017	10.35.003	10.35.001	10.35.002	10.35.004	10.35.018	10.35.007	10.35.007	Replacement.022" (10 UN)



2.9mm	2.9mm	3.2mm	2.3mm	2.3mm	2.3mm	2.3mm	3.2mm	2.9mm	2.9mm	Length
-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°	Torque
+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°	Angle
10.35.015	10.35.013	10.35.019	10.35.111	10.35.111	10.35.111	10.35.111	10.35.020	10.35.014	10.35.016	Replacement.022" (10 UN)



1 CASE • KIT 20 PIECES

SLOT22 REF 10.35.902

Prescription EDGEWISE

Morelli Slim Prescription Edgewise with hooks on canines

Slim



Length	2.5mm	2.5mm	2.5mm	2.2mm	3.0mm	3.0mm	2.2mm	2.5mm	2.5mm	2.5mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement .022" (10 UN)	10.65.206	10.65.206	10.65.204	10.65.203	10.65.201	10.65.201	10.65.203	10.65.205	10.65.206	10.65.206
Replacement .018" (10 UN)	10.65.106	10.65.106	10.65.104	10.65.103	10.65.101	10.65.101	10.65.103	10.65.105	10.65.106	10.65.106

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

Length	2.5mm	2.5mm	2.5mm	2.2mm	2.2mm	2.2mm	2.2mm	2.5mm	2.5mm	2.5mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement .022" (10 UN)	10.65.206	10.65.206	10.65.205	10.65.203	10.65.203	10.65.203	10.65.203	10.65.204	10.65.206	10.65.206
Replacement .018" (10 UN)	10.65.106	10.65.106	10.65.105	10.65.103	10.65.103	10.65.103	10.65.103	10.65.105	10.65.106	10.65.106

1 CASE • KIT 20 PIECES	
SLOT22	REF 10.65.901
SLOT18	REF 10.65.900

Morelli Slim Prescription Edgewise

Slim



Length	2.5mm	2.5mm	2.5mm	2.2mm	3.0mm	3.0mm	2.2mm	2.5mm	2.5mm	2.5mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement .022" (10 UN)	10.65.206	10.65.206	10.65.206	10.65.203	10.65.201	10.65.201	10.65.203	10.65.206	10.65.206	10.65.206
Replacement .018" (10 UN)	10.65.106	10.65.106	10.65.106	10.65.103	10.65.101	10.65.101	10.65.103	10.65.106	10.65.106	10.65.106

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

Length	2.5mm	2.5mm	2.5mm	2.2mm	2.2mm	2.2mm	2.2mm	2.5mm	2.5mm	2.5mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement .022" (10 UN)	10.65.206	10.65.206	10.65.206	10.65.203	10.65.203	10.65.203	10.65.203	10.65.206	10.65.206	10.65.206
Replacement .018" (10 UN)	10.65.106	10.65.106	10.65.106	10.65.103	10.65.103	10.65.103	10.65.103	10.65.106	10.65.106	10.65.106

1 CASE KIT 20 PIECES	
SLOT22	REF 10.65.903
SLOT18	REF 10.65.902

Morelli Standard Prescrição Edgewise

Slim



Length	3.0mm	3.0mm	3.0mm	2.3mm	3.4mm	3.4mm	2.3mm	3.0mm	3.0mm	3.0mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Ângulo	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement .022" (10 UN)	10.30.208	10.30.208	10.30.208	10.30.209	10.30.201	10.30.201	10.30.209	10.30.208	10.30.208	10.30.208
Replacement .018" (10 UN)	10.30.106	10.30.106	10.30.106	10.30.107	10.30.101	10.30.101	10.30.107	10.30.106	10.30.106	10.30.106


15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

Length	3.0mm	3.0mm	3.0mm	2.3mm	3.4mm	3.4mm	2.3mm	3.0mm	3.0mm	3.0mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Ângulo	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement .022" (10 UN)	10.30.208	10.30.208	10.30.208	10.30.209	10.30.209	10.30.209	10.30.209	10.30.208	10.30.208	10.30.208
Replacement .018" (10 UN)	10.30.106	10.30.106	10.30.106	10.30.107	10.30.107	10.30.107	10.30.107	10.30.106	10.30.106	10.30.106

1 CASE • KIT 20 PIECES	
SLOT22	REF 10.30.901
SLOT18	REF 10.30.900

Morelli Slim Prescription Ricketts with hooks on canines and premolars

2.5mm	2.5mm	2.5mm	2.2mm	3.0mm	3.0mm	2.2mm	2.5mm	2.5mm	2.5mm	Length
0°	0°	+7°	+14°	+22°	+22°	+14°	+7°	0°	0°	Torque
0°	0°	+5°	+10°	+3°	+3°	+10°	+5°	0°	0°	Angle
10.65.104	10.65.104	10.25.305	10.25.303	10.25.301	10.25.302	10.25.304	10.25.306	10.65.105	10.65.105	Replacement .018" (10 un)



2.5mm	2.5mm	2.5mm	2.2mm	2.2mm	2.2mm	2.2mm	2.5mm	2.5mm	2.5mm	Length
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Torque
0°	0°	+5°	0°	0°	0°	0°	+5°	0°	0°	Angle
10.65.105	10.65.105	10.25.307	10.65.103	10.65.103	10.65.103	10.65.103	10.25.308	10.65.104	10.65.104	Replacement .018" (10 un)

Slim




1 CASE • KIT 20 PIECES

SLOT18 REF 10.25.905

Morelli Actual Prescription Ricketts with hooks on canines and premolars

2.5mm	2.5mm	2.5mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm	Length
0°	0°	+7°	+14°	+22°	+22°	+14°	+7°	0°	0°	Torque
0°	0°	+8°	+10°	+3°	+3°	+10°	+8°	0°	0°	Angle
10.25.109	10.25.107	10.25.105	10.25.103	10.25.101	10.25.102	10.25.104	10.25.106	10.25.108	10.25.110	Replacement .018" (10 un)



2.5mm	2.5mm	2.7mm	2.2mm	2.2mm	2.2mm	2.2mm	2.7mm	2.5mm	2.5mm	Length
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Torque
0°	0°	+2°	0°	0°	0°	0°	+2°	0°	0°	Angle
10.25.110	10.25.115	10.25.113	10.25.111	10.25.111	10.25.111	10.25.111	10.25.114	10.25.116	10.25.109	Replacement .018" (10 un)

Actual




1 CASE • KIT 20 PIECES

SLOT18 REF 10.25.901

Morelli Actual Special Prescription Ricketts with hooks on canines and premolars

2.5mm	2.5mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm	Length
0°	0°	+7°	+10°	+22°	+22°	+10°	+7°	0°	0°	Torque
0°	0°	+8°	+10°	+3°	+3°	+10°	+8°	0°	0°	Angle
10.25.125	10.25.107	10.25.119	10.25.123	10.25.101	10.25.102	10.25.124	10.25.120	10.25.108	10.25.126	Replacement .018" (10 un)




2.5mm	2.5mm	2.7mm	2.2mm	2.2mm	2.2mm	2.2mm	2.7mm	2.5mm	2.5mm	Length
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	Torque
0°	0°	+2°	0°	0°	0°	0°	+2°	0°	0°	Angle
10.25.115	10.25.115	10.25.121	10.25.111	10.25.111	10.25.111	10.25.111	10.25.122	10.25.116	10.25.116	Replacement .018" (10 un)

Actual Especial



1 CASE • KIT 20 PIECES

SLOT18 REF 10.25.900

Brackets	Prescription	Use	Length.	Torque	Qty	Bonding REF SLOT 22	Bonding REF SLOT 18	Welding REF SLOT 22	Welding REF SLOT 18
	**Edgewise/ Rickets**	12 11 | 21 22	3,0 mm	0°	10	10.30.202	10.30.102	-	-
	**Edgewise**	42 41 | 31 32	1,8 mm	0°	10	10.30.203	10.30.103	-	-
	**Edgewise/ Rickets**	12 | 22 42 41 | 31 32	2,3 mm	0°	10	10.30.205	10.30.104	-	-
	**Edgewise**	15 14 13 | 23 24 25 45 44 43 | 33 34 35	2,5 mm	0°	10	10.30.206	-	-	-
	**Edgewise/ Rickets**	15 14 13 | 23 24 25 45 44 43 | 33 34 35	2,3 mm	0°	10	10.30.207	10.30.105	-	-
	**Edgewise/ Rickets**	12* | 22* 42 41** | 31 32** *Edgewise **Edgewise / Ricketts	2,3 mm	0°	10	10.30.209	10.30.107	-	-
	**Edgewise**	16 | 26 46 | 36	3,4 mm	0°	10	10.30.210	10.30.108	-	-
	**Edgewise/ Rickets**	15 14 13 | 33 34 35 23 24 25 | 45 44 43	3,0 mm	0°	10	10.30.211 10.30.212	10.30.111 10.30.112	-	-
	**Canal Vertical**	-	1,8 mm	0°	10	10.40.201	-	-	-
	**Edgewise**	11 | 21	3,4 mm	0°	10	-	-	10.31.201	-
	**Edgewise**	15 14 13 | 23 24 25 45 44 43 | 33 34 35	2,6 mm	0°	10	-	-	10.31.206	-
	**Edgewise**	15 14 13 | 23 24 25 45 44 43 | 33 34 35	2,3 mm	0°	10	-	-	10.31.207	-

Tubes



super kit TUBES

200

TUBES FOR SIMPLE BONDING

- Slot 22
- All Quadrants
- Prescription Roth with hook



4 different tubes

20.11.211	20.11.212	20.11.221	20.11.222
Roth Tube Simples .022"	Roth Tube Simples .022"	Roth Tube Simples .022"	Roth Tube Simples .022"



Simple Tube Kit with hook Prescription Roth

Prescription	Position	Length.	Tube	Torque	Offset	Qty	Bonding REF SLOT 22
Roth	17 16 26 27 47 46 36 37	3,5 mm	-	-14° -15°	14° 4°	200	20.11.921

200

TUBES FOR SIMPLE BONDING

- Slot 22
- All Quadrants
- Prescription Edgewise



Simple Tube Kit Prescription Edgewise

Prescription	Position	Length.	Tube	Torque	Offset	Qty	Bonding REF SLOT 22
Edgewise	17 16 26 27 47 46 36 37	5.2 mm	-	0°	0°	200	20.31.920

The Morelli self-ligating braces proven quality, **now in tube shape.**

Far beyond the self-ligating technique, Morelli Self-Ligating Tube combines the ease of inserting the wire into convertible tubes with the low-friction mechanics of its Passive Clip.



SLTube



SLTUBE MAIN FEATURES

- Shows Torque and Distal Offset at the base of the piece;
- Big foot base with micro-pins to ensure excellent retention during treatment;
- Easy to insert the wires, as there is no need to fit the wire into the slot;
- **Versatility:** can be used on prescriptions Capelozza, MTB and Roth;
- **Nitinol clip:** resistance and elastic memory;
- **Practicality:** no need for applying ligature;
- **Economy:** save chair time by opening and closing the clip.

Replacements on page 36

KIT

Description
SLTube Case .022" Bond - 4 quadrants
Position
17 16 | 26 27 47 46 | 36 37
REF
20.13.900



Prescription ◉ ROTH

Simple self-ligating tube



Prescription	Position	Length.	Tube	Torque	Offset	Qty	Bonding REF SLOT 22	Bonding REF SLOT 18	Welding REF SLOT 22	Welding REF SLOT 18
*	17 16	3,5 mm	–	-10°	10°	05	20.13.211	–	–	–
*	26 27	3,5 mm	–	-10°	10°	05	20.13.212	–	–	–
*	36 37	3,5 mm	–	-25°	0°	05	20.13.221	–	–	–
*	47 46	3,5 mm	–	-25°	0°	05	20.13.222	–	–	–

*Compatible with Roth, M.B.T. and Capelozza prescriptions.

Simple tube with hook



Roth	17 16	3,5 mm	–	-14°	14°	10	20.11.211	–	–	–
Roth	26 27	3,5 mm	–	-14°	14°	10	20.11.212	–	–	–
Roth	36 37	3,5 mm	–	-30°	4°	10	20.11.222	–	–	–
Roth	47 46	3,5 mm	–	-30°	4°	10	20.11.221	–	–	–

Simple tube with hook



Roth	17 16	3,5 mm	–	-14°	14°	10	–	20.11.111	20.10.211	–
Roth	26 27	3,5 mm	–	-14°	14°	10	–	20.11.112	20.10.212	–
Roth	36 37	3,5 mm	–	-30°	4°	10	–	20.11.122	20.10.222	–
Roth	47 46	3,5 mm	–	-30°	4°	10	–	20.11.121	20.10.221	–

Simple mini tube with hook



Roth	17	2,5 mm	–	-14°	0°	10	20.11.213	20.11.113	–	–
Roth	27	2,5 mm	–	-14°	0°	10	20.11.214	20.11.114	–	–
Roth	37	2,5 mm	–	-25°	0°	10	20.11.224	20.11.124	–	–
Roth	47	2,5 mm	–	-25°	0°	10	20.11.223	20.11.123	–	–

Simple convertible bracket tube



Roth M.B.T. Capelozza	44 34	2,5 mm	–	-17°	0°	10	20.11.225	–	–	–

Double convertible tube with hook



Roth	36	3,5 mm	–	-30°	4°	10	20.11.412	20.11.312	20.10.412	–
Roth	46	3,5 mm	–	-30°	4°	10	20.11.411	20.11.311	20.10.411	–

Prescription ◉ ROTH


Prescription	Position	Length.	Tube	Torque	Offset	Qty	Bonding REF SLOT 22	Bonding REF SLOT 18	Welding REF SLOT 22	Welding REF SLOT 18
Double convertible tube with hook for lip bumber										
Roth	36	3,5 mm	cervical 1,15 mm .045"	-30°	4°	10	–	–	20.10.414	20.10.314
Roth	46	3,5 mm	cervical 1,15 mm .045"	-30°	4°	10	–	–	20.10.413	20.10.313
Double convertible tube with hook										
Roth	16	3,5 mm	–	-14°	14°	10	20.11.417	–	20.10.417	–
Roth	26	3,5 mm	–	-14°	14°	10	20.11.418	–	20.10.418	–
Triple convertible tube with hook										
Roth	16	3,5 mm	oclusal 1,15 mm .045"	-14°	14°	10	20.11.611	–	20.10.611	20.10.511
Roth	26	3,5 mm	oclusal 1,15 mm .045"	-14°	14°	10	20.11.612	–	20.10.612	20.10.512
Double convertible tube with hook										
Especial	17 16	3,5 mm	–	0°	10°	10	20.11.419	–	–	–
Especial	26 27	3,5 mm	–	0°	10°	10	20.11.420	–	–	–

Compatible to Roth prescription

Prescription ◉ CAPELOZZA

Prescription	Position	Length.	Tube	Torque	Offset	Qty	Bonding REF SLOT 22	Bonding REF SLOT 18	Welding REF SLOT 22	Welding REF SLOT 18
Simple convertible tube with hook										
Capelozza Roth	17 16	3,5 mm	–	-10°	10°	10	20.19.200	–	20.19.100	–
Capelozza Roth	26 27	3,5 mm	–	-10°	10°	10	20.19.201	–	20.19.101	–
Capelozza Roth	36 37	3,5 mm	–	-25°	5°	10	20.19.203	–	20.19.103	–
Capelozza Roth	47 46	3,5 mm	–	-25°	5°	10	20.19.202	–	20.19.102	–
Double convertible tube with hook										
Capelozza (Padrão II)	17 16 | 26 27	3,5 mm	–	-10°	0°	10	20.19.301	–	20.19.300	–
Triple convertible tube with hook										
Capelozza (Padrão II)	16 | 26	3,5 mm	oclusal 1,15 mm .045"	-10°	0°	10	–	–	20.19.500	–

Prescription ◉ EDGEWISE

Prescription	Position	Length.	Tube	Torque	Offset	Qty	Bonding REF SLOT 22	Bonding REF SLOT 18	Welding REF SLOT 22	Welding REF SLOT 18
Simple convertible										
 Edgewise	17 16 / 47 46 — 26 27 / 36 37	5,2 mm	–	0°	0°	10	20.31.200	20.31.100	–	–
 Edgewise	17 16 / 47 46 — 26 27 / 36 37	5,2 mm	–	0°	0°	10	–	–	20.30.200	20.30.100
Special simple tube with hook										
 Edgewise	17 16 / 47 46 — 26 27 / 36 37	3,3 mm	–	0°	0°	10	20.31.201	–	20.30.201	–
Simple tube with hook										
 Edgewise	47 46 — 26 27	5,2 mm	–	0°	0°	10	20.31.230	–	20.30.230	–
Edgewise	17 16 — 36 37	5,2 mm	–	0°	0°	10	20.31.260	–	20.30.260	–
Simple tube with hook										
 Edgewise	17 16 / 46 47 — 26 27 / 36 37	3,5 mm	–	0°	0°	10	20.31.241	20.31.141	20.30.241	–
Simple convertible tube with hook										
 Edgewise	17 16 / 46 47 — 26 27 / 36 37	5,2 mm	–	0°	0°	10	20.31.280	20.31.180	20.30.280	20.30.180
Double tube with hook										
 Edgewise	17 16 / 46 47 — 26 27 / 36 37	3,5 mm	–	0°	0°	10	20.41.263	–	20.40.263	–
Edgewise	17 16 / 46 47 — 26 27 / 36 37	3,5 mm	–	0°	0°	10	*20.41.363		*20.40.363	

Slot .018 for occlusal and .022" for cervical
* Slot .018* for occlusal and .022" for cervical

Prescription ◉ RICKETTS

Prescription	Position	Length.	Tube	Torque	Offset	Qty	Bonding REF SLOT 22	Bonding REF SLOT 18	Welding REF SLOT 22	Welding REF SLOT 18
Double tube with hook										
 Ricketts	17 16 / 46 47 — 26 27 / 36 37	3,5 mm	–	0°	0°	10	–	20.41.163	–	20.40.163

Prescription ◉ EDGEWISE/RICKETTS

Prescription	Position	Length.	Tube	Torque	Offset	Qty	Bonding REF SLOT 22	Bonding REF SLOT 18	Welding REF SLOT 22	Welding REF SLOT 18
Double convertible tube with hook										
Edgewise Ricketts	16 26 / 46 36	3,5 mm	cervical 1,15 mm .045"	0°	0°	10	-	-	20.40.331	20.40.330
Edgewise Ricketts	16 26 / 46 36	3,5 mm	cervical 1,30 mm .051"	0°	0°	10	-	-	20.40.431	20.40.430



Double convertible tube with hook										
Edgewise Ricketts	16 26 / 46 36	3,5 mm	-	0°	0°	10	20.41.501	20.41.500	20.40.501	20.40.500



Triple tube with hook with hook										
Edgewise Ricketts	16 26 / 46 36	3,5 mm	central 1,15 mm .045"	0°	0°	10	-	-	20.50.181	20.50.180



Prescription ◉ M.B.T.

Simple convertible tube with hook										
M.B.T.	37 47	3,5 mm	-	-10°	0°	10	20.35.215	-	-	-



Simple convertible tube										
M.B.T.	17 16	3,5 mm	-	-14°	10°	10	20.35.113	-	20.35.103	-
M.B.T.	26 27	3,5 mm	-	-14°	10°	10	20.35.114	-	20.35.104	-



M.B.T.	46 36	3,5 mm	-	-20°	0°	10	20.35.212	-	20.35.202	-



Double convertible tube with hook										
M.B.T.	46 36	3,5 mm	-	-20°	0°	10	-	-	20.35.401	-



Prescription ◉ M.B.T.

Prescription	Position	Length.	Tube	Torque	Offset	Qty	Bonding REF SLOT 22	Bonding REF SLOT 18	Welding REF SLOT 22	Welding REF SLOT 18
Double convertible tube with hook										
M.B.T.	16	3,5 mm	cervical 1,15 mm .045"	-14°	10°	10	–	–	20.35.402	–
M.B.T.	26	3,5 mm	cervical 1,15 mm .045"	-14°	10°	10	–	–	20.35.403	–
Double convertible tube with hook for P.L.A										
M.B.T.	46 | 36	3,5 mm	cervical 1,15 mm .045"	-20°	0°	10	–	–	20.35.405	–
Double convertible tube with hook										
M.B.T.	46 | 36	3,5 mm	–	-20°	0°	10	20.35.416	–	–	–
Triple convertible tube with hook										
M.B.T.	16	3,5 mm	oclusal 1,15 mm .045"	-14°	10°	10	–	–	20.35.501	–
M.B.T.	26	3,5 mm	oclusal 1,15 mm .045"	-14°	10°	10	–	–	20.35.502	–

Lingual ◉ LINGUAL

Prescription	Position	Length.	Torque	Offset	Qty	Welding 2 x 0,90mm 2 x .036"
Lingual tube with hook						
Lingual	16 26 / 46 36	3,5 mm	0°	0°	10	20.61.206

Adhesives



ORTHOBOND® PLUS
Color Change Adhesive

Made by those who best understand about brackets, for those who best understand about orthodontics.

High performance

ORTHOBOND offers high performance for bonding metal, ceramic and composite brackets, tubes and other accessories.

Yields more than 15 complete cases

Color Change System

After light curing, the adhesive loses its pink tone and has a similar color to the tooth.



The use of **ORTHOBOND**PLUS together with the **ORTHOPRIMER** provides a considerable increase in **BONDING EFFICIENCY**

Yields more than 50 complete cases

ORTHOPRIMER is a fluid adhesive used as a base for activating photopolymerizing orthodontic adhesives.

It has high fluidity for micro-penetration into enamel prisms, enhancing bonding to brackets and accessories.

ORTHOPRIMER

Adhesive Orthobond Plus

Orthobond Plus Color Change Adhesive

Size	Qty	REF
4g	1	85.01.017



Cap it immediately after each use: it hardens in ambient light!



Adhesive OrthoPrimer

OrthoPrimer

Size	Qty	REF
5g	1	85.01.016





Check out essential bonding tips in our virtual Orthodontic Bonding Manual.

OPENBITE

Bite-Raising Resin

Faster disocclusion and safer removal



Secure identification

Formula with permanent blue pigment, which allows easy visualization from application to removal of occlusal raising, **avoiding unnecessary wear on tooth enamel.**

Resistance

Its high adhesion and durability provide excellent mechanical and adhesive resistance, preventing wear during chewing.

Photopolymerizable

Its high adhesion and durability provide excellent mechanical and adhesive resistance, preventing wear during chewing.

Self-adhesive

Its formula has good adhesion and may eliminate the need for a primer in some cases, speeding up the clinical stage.

Viscosity

Open Bite has the ideal consistency, which makes it versatile for building bite lifts.

Contents	Qty.	REF
4 g	1	85.01.018


Accessories



Long Stem Hook

for conventional and self-ligating metal brackets

The Long Stem Hook for Brackets is an orthodontic solution designed to optimize distal canine traction.





Made from high-strength stainless steel, the accessory offers **precision, safety and compatibility** with both self-ligating and conventional techniques. All of this with the Morelli quality standard.

Why choose Long Stem Sook for brackets?

- Direct attachment to the bracket, providing an efficient anchorage point;
- Facilitates the application of orthodontic forces at the tooth's center of resistance, optimizing the movement of the canine;
- It contributes to more predictable and controlled treatments;
- Durability and resistance thanks to stainless steel construction.

Cada código possui encaixe perfeito nos bráquetes para os quais foram desenvolvidos.

Long Stem Hook for Metal Brackets Roth Max Canines	**Long Stem Hook for Metal Brackets** Roth SLI, SLP e SLH Canines
	
REF 30.40.050	REF 30.40.051
4 units	4 units
Position: 13 | 23 / 43 | 33	

Super Kit
Orthodontic Accessories

32 **TYPES OF PRODUCT**
KIT WITH 240 PIECES






COMES WITH ORGANIZER BOX with numbered marking of all accessories


Orthodontic Accessories Kit	Type	Quantity	REF
	32 types of orthodontic accessories	240	30.90.900

Tubo Bráquete Prescrição Roth Simples Conv. 20.11.225 10 un. 01	Tubo Lingual 2x.036" Solda 16,26,36,46 20.61.206 10 un. 02	Contenção Lingual TwistForm 30.06.303 06 un. 03	Botão para Colagem Plano 30.10.205 30 un. 04	Botão Ceramic 30.10.210 05 un. 05	Gancho Soldagem 30.20.101 10 un. 06	Gancho para Tração de Dentes Inclusos 30.21.104 10 un. 07	Gancho para Tração de Dentes Inclusos com Corrente 30.21.107 05 un. 08
Gancho Ponta Bola Direito 30.40.003 10 un. 09	Gancho Ponta Bola Esquerdo 30.40.004 10 un. 10	Gancho Ponta Bola Reto 30.40.005 10 un. 11	Stop Aberto 30.40.008 10 un. 12	Gancho Ponta Bola Universal 30.40.009 10 un. 13	Gancho Universal para MPO 30.40.010 10 un. 14	Tubo Cruzado com Recuo Cervical .022" 30.40.018 10 un. 15	Stop Fechado com Aplicador 30.40.021 20 un. 16
Cursor Curto Sup. Esquerdo / Inf. Direito 30.45.011 02 un. 17	Cursor Curto Sup. Direito / Inf. Esquerdo 30.45.012 02 un. 18	Cursor Médio Sup. Esquerdo / Inf. Direito 30.45.013 02 un. 19	Cursor Médio Sup. Direito / Inf. Esquerdo 30.45.014 02 un. 20	Cursor Longo Sup. Esquerdo / Inf. Direito 30.45.015 02 un. 21	Cursor Longo Sup. Direito / Inf. Esquerdo 30.45.016 02 un. 22	Stop e Guia de Mordida 30.60.001 10 un. 23	Mini Esporão Educador Lingual 30.60.002 10 un. 24
Mola Fechada NiTi - 7mm 35.20.064 04 un. 25	Mola Fechada NiTi - 9mm 35.20.065 04 un. 26	Mola Fechada NiTi - 12mm 35.20.066 04 un. 27	Mola Fechada NiTi - 15mm 35.20.067 04 un. 28	Mola Fechada NiTi Miniparafuso 7mm 37.30.001 04 un. 29	Mola Fechada NiTi Miniparafuso 9mm 37.30.002 04 un. 30	Mola Fechada NiTi Miniparafuso 12mm 37.30.003 04 un. 31	Mola Fechada NiTi Miniparafuso 15mm 37.30.004 04 un. 32

Number of pieces per cavity in this package



Universal hook for MPO

Type	Size	Qty	REF
slot .022”		05	30.40.010

Height-adjustable hooks, remove the hooks by flexing to the desired size and apply finishing



Hook Anchoring Ancoragem

Used in the orthodontic arch to anchor elastic chains, springs, metal ligature wires, to aid treatments involving retraction, mesialization and anchorage.

They anchorage is better in rectangular arches, as they are larger in gauge.




Ball curved tip hook

Type	Size	Qty	REF
Direito Right **slot .022”**		10	30.40.003
Esquerdo Left **slot .022”**			30.40.004




Ball straight tip hook

Type	Size	Qty	REF
slot .022”		10	30.40.005




Ball universal tip hook

Type	Size	Qty	REF
slot .022”		10	30.40.009



Hook for traction of non erupted teeth with chain


Size	Qty	REF
	10	30.21.107

Hook for traction


The hook is bonded to the impacted tooth and the chain remains anchored in the intraoral arch. The tooth extrusion occurs as the links appear, which, once counted, identify the movement and positioning of the tooth.




Hook for welding

Size	Qty	REF
	100	30.20.101

Hook for traction of non erupted teeth

Size	Qty	REF
	10	30.21.104



Gurin

These are devices applied to orthodontic arches with the aid of a hexagonal key*, allowing devices to be anchored directly to the arch.

*Gurin key page 114 - 75.02.011


Gurin with closed round slot

Type	Size	Qty	REF
Universal	.022” adjustable	10	35.05.001



Gurin with open rectangular round hook

Type	Size	Qty	REF
Universal	.022” adjustable	10	35.05.011



Gurin with closed round slot and hook

Type	Size	Qty	REF
Right	.022” adjustable	10	35.05.002
Left			35.05.003



Gurin with open rectangular slot and hook

Tipo *Type*	Medidas *Size*	Qtd *Qty*	REF
Right	.022” adjustable	10	35.05.012
Left			35.05.013



How to use gurin o Gurin

Gurin with closed round slot and right hook

Type	Size	Qty	REF
Right	.022” adjustable	05	35.05.004
Left			35.05.005



Button

Accessory used in tooth movement such as: tooth rotation, traction, retraction or anchorage systems.
They can trap elastics and ligature wires.

When to use the orthodontic button

Button Welding



Type	Size	Qty	REF
For Welding		10	30.10.105


Button for Bonding Convex



Type	Size	Qty	REF
Convex for bonding colagem		10	30.10.202


Button for Bonding Concave



Type	Size	Qty	REF
Convex for bonding colagem		10	30.10.203


Button for Bonding Flat



Type	Size	Qty	REF
For bonding flat		10	30.10.205


Ceramic button



Type	Size	Qty	REF
Ceramic for bonding		02	30.10.210


PolySafira Button



Type	Size	Qty	REF
Sapphire for bonding		02	30.10.211

Composite button




Type	Size	Qty	REF
Composite for bonding		10	30.10.206

Power Arm Ceramic

Ideal for use in mesialization, distalization, space closure, and anterior retraction mechanics.
Single-body accessory, injected from translucent alumina ceramic.

Power arm ceramic for mesialization or distalization

Size	Qty	REF
	01	30.21.001




Stop

For rectangular arches: Efficient for containing or limiting the movement of arches and orthodontic accessories, as well as avoiding the discomfort caused by the displacement of the arch at the distal end of the tube.


For round arches: Designed for demarcating the midline or locking round orthodontic arch, mainly to avoid discomfort due to displacement of the arch at the distal end of the tube. Limits lateral movement without causing deformation. It is used in all orthodontic techniques, including self-ligating brackets.


Stop for rectangular arches

Size	Qty	REF
 .022”	10	30.40.008




Stop for round arches

Size	Qty	REF
 .021” / Ø 0.54mm	25	30.40.020



Stop closed with applicator for wires up to .022"

Size	Qty	REF
 .021” / Ø 0.54mm	10	30.40.021




Check how to apply aplicar

Stop and Guide

The Stop and Guide main function is to reduce the impact of the upper teeth on the lower brackets, as it works by opening up the bite, while at the same time re-educating and resting the muscles.

Stop and Bite Guide


Size	Qty	REF
 Base / Base **3.6mm x 4.0mm**	10	30.60.001



Mini Spur

The mini spur induces the tongue to reduce pressure on the anterior, upper and/or lower teeth, helping to treat digital sucking (thumb sucking), correct open bites and re-educate relapses.

Mini spur tongue educator

Size	Qty	REF
 Base / Base **2.8mm x 3.8mm**	10	30.60.002



Spring

Morelli springs offer light and continuous forces, allowing you to increase the activation intervals. They are excellent devices for obtaining tensile or compressive forces. The closed springs for traction have ends prepared for adapting to accessories.



Open Spring CrNi

Type	Length	Size	Qty	REF
Open compression	1 Meter	 ø interno / internal 0,80mm (.032")	01	35.20.060




Closed Spring CrNi

Type	Length	Size	Qty	REF
Closed traction	1 Meter	 ø interno / internal 0,80mm (.032")	01	35.20.061



Spring NiTiS

Type	Length	Size	Qty	REF
Open compression	178 mm	 ø interno / internal 0,76mm (.030")	03	35.20.062



Tip: To avoid the open spiral tip, it is recommended to cut the spring at the closed spiral.

Spring NiTi open and closed


Type	Length	Size	Qty	REF
Open and closed compression	178 mm	 ø interno / internal 1,00mm (.040")	03	35.20.072



ø 1.6 mm

Spring NiTi closed

Type	Length	Range of application	Size	Qty	REF
Closed traction	7 mm	10 a 15 mm	200gf	04	35.20.064
	9 mm	15 a 25 mm			35.20.065
	12 mm	25 a 35 mm			35.20.066
	15 mm	35 a 45 mm			35.20.067



ø 2.3 mm

Spring NiTi closed for mini screw

Type	Length	Range of application	Size	Qtd *Qty*	REF
Closed traction	7 mm	10 a 15 mm	200gf	04	37.30.001
	9 mm	15 a 25 mm			37.30.002
	12 mm	25 a 35 mm			37.30.003
	15 mm	35 a 45 mm			37.30.004

Cursor

Morelli Sliding Jig is an intraoral device for installation in the intraoral arch, which allows springs or elastics to be anchored to an adjustable position hook. The clinical force is adjusted by positioning the hook along the Ø 0.7 mm wire sliding jig, using the Gurin* key.

**Gurin key page 114

Sliding jig for Mesialization and Distalization

Type	Size	Qty	REF
12mm Ø1,30mm	Short - Sup. Esquerdo Inf. Direito Short - UL/LR	05	30.45.011
12mm	Short - Sup. Esquerdo Inf. Esquerdo Short - UR/LL		30.45.012
18mm	Medium- Sup. Esquerdo Inf. Direito Medium - UL/LR	05	30.45.013
18mm	Medium- Sup. Direito Inf. Esquerdo Medium - UR/LL		30.45.014
32mm	Long- Sup. Esquerdo Inf. Direito Long - UL/LR	05	30.45.015
32mm	Long- Sup. Direito Inf. Esquerdo Long - UR/LL		30.45.016


Fio Ø 0,7mm








Distalization



Mesialization



Check out the application of the sliding jigs

Retention

Orthodontic retention is used to keep teeth in their new positions after orthodontic movement with the substantial function of providing the stability of the results obtained after the end of orthodontic treatment.

TwistForm Retainer 3x3

Type	Length	Size	Qty	REF
3 wires	41,5 mm	Ø externo / esternal 0,50mm (.020")	10	30.06.303
3 wires	41,5 mm	Ø externo / esternal 0,65mm (.026")	10	30.06.304



Check out the TwistForm application



Tube for activators bar

Size	Qty	REF
Ø 1,15mm Ø 0.45”	10	30.05.145
Ø 1,30mm Ø 0.51”	10	30.05.151



Cross Tube

Tipy	Size	Qty	REF
-	.022”	10	30.40.006

Cross Tube Cruzado

The cross tube is attached to the orthodontic arch to receive wires up to .021"x.025". Coupled to the main arch, it allows the installation of segmented arches, cantilevers, straps and others. It helps with space closure, anterior segment intrusion and loss of anchorage.

Cross tube with cervical offset: The tube cervical offset prevents the antagonist tooth from touching the wire inserted into the hole.



Cervical Offset cross tube

Tipy	Size	Qty	REF
-	.022”	10	30.40.016

”

Molar distalizer type Jones universal

Size (L)	Compression*	Force	Qty	REF
15 mm	3 mm	60gf	02	30.45.003
	4 mm	65gf		
	5 mm	70gf		
	6 mm	75gf		
	7 mm	80gf		


Distalizer

This device is used to distalize upper molars and offers constant forces. Activation of the nitinol spring releases forces of 60 to 80 gf when subjected to average compressions of 5 mm.

Mola NiTi



Movable cervical/ occlusal eyelet for anchorage




*Medium Values

Telescopic Tube

Made of stainless steel, they are used in the composition of some orthopaedic and orthodontic appliances that feature sliding movements on a wire or stem. For example, Mandibular Protractor Appliances (MPAs) are used to correct class II malocclusions

Telescopic tube

Indicated for wires up to		Qty	Length	REF
Internal diameter	0,80mm	05	330mm	30.40.030
Internal diameter	0,90mm	05	330mm	30.40.031
Internal diameter	1,40mm	05	330mm	30.40.032



Headed Ball Clasp

The Headed Ball Clasp is used in the manufacture of removable appliances to create "retention" points by means of its ball, which presses against the teeth.

Headed ball clasp

Dimensions	Qty	Length	REF
Stem: Ø 0,8mm Ball: Ø 1,5mm	50	38mm	30.40.041





Palatal bar

The palatal bar connects to the lingual tubes. Used for a variety of purposes, it can generate forces for tooth rotation, intrusion or torque, as well as acting as anchor reinforcement.

Palatal bar with mesial and distal loop

Size		Qty	REF - Mesial	REF - Distal
*34mm	Ø 0,9mm	10	30.30.034	30.30.134
37mm			30.30.037	30.30.137
39mm			30.30.039	30.30.139
41mm			30.30.041	30.30.141
43mm			30.30.043	30.30.143
45mm			30.30.045	30.30.145
47mm			30.30.047	30.30.147
49mm			30.30.049	30.30.149
52mm			30.30.052	30.30.152
55mm			30.30.055	30.30.155





Mire mesh 80

Type	Size	Qty	REF
CrNi	Fio = Ø 0,12mm 100 x 100 mm	01	35.11.001

Wire Mesh

Used to make the base for bonding tubes and other pieces. Uses 80 stainless steel stranded with 80 wires/inch





Silver Solder

Type	Size	Qty	REF
56% (Ag) Fusion = 650°C	Roll Ø 0,05mm 5 m	01	35.30.010
	Stick Ø 0,8mm 300 m	04	35.30.011
	Stick Ø 0,8mm 300 m	10	35.30.013

Silver Solder

Easy-to-apply silver welder free of toxic materials.



Solder flux

Size	Qty	REF
20 ml	01	35.30.020

Pasta para solda

Material utilizado como fluxo indispensável para unir metais através de brasagem, permitindo soldas uniformes e sem bolhas.



Orthodontic wax protection

Qty	REF
10 cards	35.33.001

Cera Ortodôntica

A Cera de Proteção desenvolvida pela Morelli oferece conforto durante a adaptação ao aparelho ortodôntico, além de evitar as consultas de urgência e aumentar a colaboração dos pacientes em início do tratamento. Aromatizada com suave sabor de menta.



Orthodontic Wax

Qty	REF
1 cards	35.33.002

Aesthetic Solutions



Technology and High Performance for treatments that require discretion

The **Polysafira Button** and **Power Arm Ceramic** are solutions that complement each other to provide greater aesthetics, adhesion and movement control in treatments with aesthetic brackets or aligners.



Polysafira Button
REF 30.10.211
2 units

Power Arm Ceramic
REF 30.21.001
1 unit

Aesthetic Discretion: Made from translucent alumina ceramic

Movement Control: It allows forces to be applied at two different heights, enabling precise mesialization or distalization movements with aligners or brackets.

Monobloc Design: Injection-molded, single-body accessory that guarantees strenght and practicality when bonding, ideal for canines and first premolars.

Integration with Anchoring Systems: It transmits forces effectively when used with mini screws and other systems, offering biomechanical versatility.

Upper Aesthetics: Transparency that adapts to the natural color of the teeth.

Excellent Adhesion: Bi-blasted, marked base, with mechanical and chemical retention compatible with the main orthodontic adhesives on the market.

Injection Technology: Produced in polycrystalline ceramic using the C.I.M. (Ceramic Injection Molding) process, which guarantees precision and resistance.

Versatility: Ideal for treatment with aesthetic brackets and aligners.

Intelligent Design: Optimized anatomical base for better adaptation in more restricted areas, such as the dental crowns cervical regions.

CLARIS
THERMOPLASTIC SOLUTIONS



Your knowledge in control of new possibilities

Claris Plates have been developed for thermoforming in numerous dental applications.

They are resistant and safe for malleability by vacuum or pressure heating. In addition to excellent final shape memory, with high ductility, it does not fracture under the forces applied during use.

Dental applications:

- Aesthetic Aligners
- Retention plates
- Indirect bonding guides
- Occlusal plates
- Trays for whitening or fluoride application



AVAILABLE SIZES

Size	Qty.	REF
Plate 0,50mm - Ø 125mm	10	66.01.050
Plate 0,76mm - Ø 125mm	10	66.01.076
Plate 1,00mm - Ø 125mm	10	66.01.100

Check out how to use Claris Plates

See other **complimentary products**



Turquoise Blue
80.04.216

Pocket-size Boxes

As well as being beautiful, Morelli Pocket-size Boxes are super resistant, practical and hygienic.

They have a vent at the top of the lid, which ensures air circulation, keeping both the box and its contents dry at all times.

Suitable for storing and transporting retaining, occlusal, whitening plates, mobile appliances and aligners.



Navy Blue
80.04.217

Available in 10 colors:



Blue / White
80.04.211

Pink / White
80.04.212

Off White
80.04.214

Piano Black
80.04.215

Cherry Red
80.04.218

Magenta Pink
80.04.219

Tangerine Orange
80.04.220

Lime Green
80.04.221

Orthodontic positioner and extractor

This instrument facilitates the placement and removal of elastics (intraoral and extraoral).

It also helps to remove aligners and retainers safely, as it avoids direct contact with the nails, preventing possible breakage and contamination.



Code **75.01.015**



Code **30.10.210**

Ceramic button


A versatile piece with superior aesthetics and resistance. Allows for special movements, especially in aligner treatments.

Recently the design has been improved, allowing bonding in narrow areas such as the cervicals of teeth.

Absolute Anchoring



A complete range of **orthodontic anchorage**



Main features of Morelli mini screws

- **Active tip with self-drilling profile,** specially designed to offer stability, a secure fit in the bone texture and resistance for immediate orthodontic load applications;
- **Transmucosal section with perfect adaptation** and maintenance of the peri-implant soft tissues;
- **High biocompatibility** material tested and approved by ANVISA for biological safety;
- **Ø 0.8 mm hole (suitable for wires up to .021”x.025”)** in the hexagonal head, specially dimensioned for the orthodontic accessories application;
- **Easy to use:** all Mini screws are compatible with the same Morelli Tool Kit;
- **Availability:** available from all major dental distributors in Brazil and guaranteed to be available at any time;
- **Fair prices:** in order to contribute to the evolution of Brazilian orthodontics, we try to offer affordable prices.



Titanium Mini Screw | Button Head

Ideal for interradicular application, it allows for various possibilities of movement during orthodontic treatment, regardless of the patient’s cooperation.



Titanium Mini Screw | Bracket Head

In addition to the spring or ligatures anchoring, it allows the use of cantilevers, and, therefor, it offers new possibilities in orthodontic movement. The Bracket Head Mini screw has two 0.021” x 0.025” slots perpendicular to each other in the screw head to precisely accommodate rectangular wires.



Steel Mini Screw | Button Head


Ideal for extra-radicular application, the Steel Mini screw, when compared to titanium, offers greater resistance to bending and torsion forces.




Mini Screw with Rectangular Hole | Button Head

Expanding the possibilities for tooth movement and offering even more solutions for correcting malocclusions, the Rectangular Hole Orthodontic Mini screws have been developed in two models: steel, for extra-radicular installation, and titanium, for interradicular installation


Titanium Orthodontic Mini screw with Button Head



Threading Length	Transmucosal	Diameter	Insertion	Qty	REF
6mm	1mm	Ø1.5mm	Clockwise	2	37.10.101
8mm					37.10.102
6mm	2mm	Ø1.5mm	Clockwise	2	37.10.201
8mm					37.10.202
10mm					37.10.203




Titanium Orthodontic Mini screw with Bracket Head



6mm	1mm	Ø1.5mm	Counter	1	37.10.231
8mm					37.10.232
6mm	1mm	Ø1.5mm	Clockwise	1	37.10.241
8mm					37.10.242




Titanium Orthodontic Mini screw Rectangular Hole




8mm	2mm	Ø1.5mm	Clockwise	1	**37.10.210**


Steel Orthodontic Mini screw Rectangular Hole



10mm	6mm	Ø2.0mm	Clockwise	1	**37.10.603**



Steel Orthodontic Mini screw with Button Head



5mm	4mm	Ø2.0mm	Clockwise	1	37.10.405
7mm	4mm	Ø2.0mm	Clockwise	1	37.10.407
10mm	4mm	Ø2.0mm	Clockwise	1	37.10.410
12mm	4mm	Ø2.0mm	Clockwise	1	37.10.412
14mm	4mm	Ø2.0mm	Clockwise	1	37.10.414



All Morelli Mini screws are sterilized.


Spear Tip

Used to facilitate entry in the gums and marking the bone to be drilled by the Mini screw.



121°C

Spear Tip

Size	Qty	REF
	1	37.20.010


Hexagonal Key

Long, medium and short to fit the mini screw. So, it's equipped with a Nitinol retention ring, durable and safe.



121°C


Hexagon Key Long Stem

Size	Qty	REF
	1	37.20.011



121°C


Hexagon Key Medium Stem

Size	Qty	REF
	1	37.20.015



121°C

Hexagonal Key Short Stem

Size	Qty	REF
	1	37.20.012

Tool Kit

The excellence and functionality you already know, now also in the **pink version.**



The Morelli Tool Kit includes:



Spear Tip
REF 37.20.010



Long Hexagonal Key
REF 37.20.011



Medium Hexagonal Key
REF 37.20.015



Short Hexagonal Key
REF 37.20.012




Digital Hexagonal Key
REF 37.20.013



Hand Key
ROSA - REF 37.20.017
AZUL - REF 37.20.014


Color	Qty.	REF
Pink	1	37.20.003
Blue	1	37.20.001



Digital hexagonal key

This screwdriver can be used in confined spaces to install and remove of orthodontic mini screws.

Hexagonal Digital Key


Size	Qty	REF
	1	37.20.013



Contra Angle Key

Coupled to the dental contra-angle, the Key replaces the low-speed micromotor, transforming the contra-angle into an angled installation key.

Contra Angle Key


SIZE	QtY	REF
	1	37.20.016




Orthodontic Drill

To increase safety during the mini screws installation, in certain cases it is recommended to pre-drill with mini screw bits, as this facilitates access to the installation site and reduces the insertion torque. Indicated as an auxiliary in the installation of mini screws, Orthodontic Drills contribute to a safe and predictable procedure, especially in situations of dense bone or difficult access.

Orthodontic Drill for Mini screw Ø 1,1 mm - Short

Size	Qty	REF
	1	37.20.018



Size	Qty	REF
	1	37.20.019



Bands



Universal band for molar

Universal bands have internal perimeters ranging from 28 to 42 mm, corresponding to the perimeters of the teeth.
.



Band for 1st and 2nd molar lower right/left

Band for 1st and 2nd molar upper right/left

Perimeter mm	Qty	REF	
		47 46 36 37	17 16 26 27
28	10	40.01.280	40.02.280
28,5	10	40.01.285	40.02.285
29	10	40.01.290	40.02.290
29,5	10	40.01.295	40.02.295
30	10	40.01.300	40.02.300
30,5	10	40.01.305	40.02.305
31	10	40.01.310	40.02.310
31,5	10	40.01.315	40.02.315
32	10	40.01.320	40.02.320
32,5	10	40.01.325	40.02.325
33	10	40.01.330	40.02.330
33,5	10	40.01.335	40.02.335
34	10	40.01.340	40.02.340
34,5	10	40.01.345	40.02.345
35	10	40.01.350	40.02.350
35,5	10	40.01.355	40.02.355
36	10	40.01.360	40.02.360
36,5	10	40.01.365	40.02.365
37	10	40.01.370	40.02.370
37,5	10	40.01.375	40.02.375
38	10	40.01.380	40.02.380
38,5	10	40.01.385	40.02.385
39	10	40.01.390	40.02.390
39,5	10	40.01.395	40.02.395
40	10	40.01.400	40.02.400
40,5	10	40.01.405	40.02.405
41	10	40.01.410	40.02.410
41,5	10	40.01.415	40.02.415
42	10	40.01.420	40.02.420



Anatomical band for molar

Morelli anatomical bands reproduce the morphology of molars, allowing for optimum adaptation. Available for upper and lower molars, right and left side.



Lower right/left

Upper right/left

Perimeter mm	Qty	REF	
		47 46 36 37	17 16 26 27
29	10	40.12.290	40.22.290
29,5	10	40.12.295	40.22.295
30	10	40.12.300	40.22.300
30,5	10	40.12.305	40.22.305
31	10	40.12.310	40.22.310
31,5	10	40.12.315	40.22.315
32	10	40.12.320	40.22.320
32,5	10	40.12.325	40.22.325
33	10	40.12.330	40.22.330
33,5	10	40.12.335	40.22.335
34	10	40.12.340	40.22.340
34,5	10	40.12.345	40.22.345
35	10	40.12.350	40.22.350
35,5	10	40.12.355	40.22.355
36	10	40.12.360	40.22.360
36,5	10	40.12.365	40.22.365
37	10	40.12.370	40.22.370
37,5	10	40.12.375	40.22.375
38	10	40.12.380	40.22.380
38,5	10	40.12.385	40.22.385
39	10	40.12.390	40.22.390
39,5	10	40.12.395	40.22.395
40	10	40.12.400	40.22.400
40,5	10	40.12.405	40.22.405
41	10	40.12.410	40.22.410
41,5	10	40.12.415	40.22.415
42	10	40.12.420	40.22.420

5 pieces RL + 5 pieces LL



5 pieces RU + 5 pieces LU



Band for upper premolar

BU

Upper

Premolar bands are anchorage devices used in maxillary expansion treatments associated with Haas and Hyrax breakers and distalizers, Jones Jig and Nance buttons. It can also be used in cases of recurrent bracket falls, where the bracket is welded to the band.

Size	Qty	REF 15 14 | 25 24
12	10	40.50.012
13	10	40.50.013
14	10	40.50.014
15	10	40.50.015
16	10	40.50.016
17	10	40.50.017
18	10	40.50.018
19	10	40.50.019
20	10	40.50.020
21	10	40.50.021
22	10	40.50.022
23	10	40.50.023
24	10	40.50.024
25	10	40.50.025
26	10	40.50.026
27	10	40.50.027
28	10	40.50.028
29	10	40.50.029
30	10	40.50.030
31	10	40.50.031
32	10	40.50.032





See bands seater on page 107

Bands with tubes for delivery immediate.

To order your Anatomical Band with Welded Tube, follow these steps:

Step1

Choose the "Lower or Upper" band size.

#TIP

Visit the page to place your order

Upper Anatomical Band:				TUBOS					
				SIMPLES	DUPLO		TRIPLO	TRIPLO + LINGUAL	
Prescrição				ROTH	ROTH	M.B.T. "L.B.	ROTH	ROTH	LINGUAL
Torque / Torque				-14°	-14°	-14°	-14°	-14°	–
Distal / Offset				14°	14°	10°	14°	14°	–
Tubo / Tube ø				–	–	.045"	.045"	.045"	–
Slot				.022"	.022"	.022"	.022"	.022"	2x .036"
Perímetro Perimeter mm	Bandas Bands	R	Qtd	20.10.211	20.10.417	20.35.402	20.10.611	20.10.611	20.61.206
		L	Qty	20.10.212	20.10.418	20.35.403	20.10.612	20.10.612	–
32	40.22.320	R L	2	47.11.007	–	–	–	–	
32,5	40.22.325	R L	2	47.11.008	–	–	–	–	
33	40.22.330	R L	2	47.11.009				48.31.109	
33,5	40.22.335	R L	2	47.11.010	–	–	47.31.010	48.31.110	
34	40.22.340	R L	2	40.11.011	47.21.011	47.21.111	47.31.011	48.31.111	
34,5	40.22.345	R L	2	47.11.012	47.21.012	47.21.112	47.31.012	48.31.112	
35	40.22.350	R L	2	47.11.013	47.21.013	47.21.113	47.31.013	48.31.113	
35,5	40.22.355	R L	2	47.11.014	47.21.014	47.21.114	47.31.014	48.31.114	
36	40.22.360	R L	2	47.11.015	47.21.015	47.21.115	47.31.015	48.31.115	
36,5	40.22.365	R L	2	47.11.016	47.21.016	47.21.116	47.31.016	48.31.116	
37	40.22.370	R L	2	47.11.017	47.21.017	47.21.117	47.31.017	48.31.117	
37,5	40.22.375	R L	2	47.11.018	47.21.018	47.21.118	47.31.018	48.31.118	
38	40.22.380	R L	2	47.11.019	47.21.019	47.21.119	47.31.019	48.31.119	
38,5	40.22.385	R L	2	47.11.020	47.21.020	47.21.120	47.31.020	48.31.120	
39	40.22.390	R L	2	47.11.021	47.21.021	47.21.121	47.31.021	48.31.121	
39,5	40.22.395	R L	2	47.11.022	–	–	47.31.022	48.31.122	
40	40.22.400	R L	2	47.11.023	–	–	47.31.023	48.31.123	
40,5	40.22.405	R L	2	47.11.024	–	–	47.31.024	48.31.124	
41	40.22.410	R L	2	47.11.025	–	–	47.31.025	48.31.125	
41,5	40.22.415	R L	2	47.11.026	–	–	47.31.026	48.31.126	
42	40.22.420	R L	2	47.11.027	–	–	47.31.027	48.31.127	

Step 2
Choose the Tube


Step 3
Identify the crossing code to place your order

Examples:

You will receive 1 pair of the chosen set (Right/Left).

Family 47
Band with Tube

Right(R) Left(L)



Band 40.12.320



Tube 20.10.221 Tube 20.10.222





Code
47.10.007
2 un

Family 48
Band with Tube

Right(R) Left(L)



Band 40.22.330



Tube 20.10.611 Tube 20.10.612 Tube 20.61.206







Code
48.31.109
2 un

ANATOMIC BAND with tube Lower				TUBES					
				SIMPLE	DOUBLE			DOUBLE+ LINGUAL	
Prescription				ROTH	ROTH	ROTH *L.B.	M.B.T.	ROTH	LINGUAL
Torque				-30°	-30°	-30°	-20°	-30°	-
Offset				4°	4°	4°	0°	4°	-
Tube ø				-	-	.045"	-	-	-
Slot				.022"	.022"	.022"	.022"	.022"	2x .036"
									
Perimeter mm	Bands	R / L	Qty	20.10.221 / 20.10.222	20.10.411 / 20.10.412	20.10.413 / 20.10.414	20.35.401 / -	20.10.411 / 20.10.412	20.61.206 / 20.61.206
32	40.12.320	R L	2	47.10.007	47.20.007	-	-	-	
32,5	40.12.325	R L	2	47.10.008	47.20.008	-	-	-	
33	40.12.330	R L	2	47.10.009	47.20.009	-	-	-	
33,5	40.12.335	R L	2	47.10.010	47.20.010	-	-	-	
34	40.12.340	R L	2	47.10.011	47.20.011	-	-	-	
34,5	40.12.345	R L	2	47.10.012	47.20.012	47.20.112	47.20.312	48.20.112	
35	40.12.350	R L	2	47.10.013	47.20.013	47.20.113	47.20.313	48.20.113	
35,5	40.12.355	R L	2	47.10.014	47.20.014	47.20.114	47.20.314	48.20.114	
36	40.12.360	R L	2	47.10.015	47.20.015	47.20.115	47.20.315	48.20.115	
36,5	40.12.365	R L	2	47.10.016	47.20.016	47.20.116	47.20.316	48.20.116	
37	40.12.370	R L	2	47.10.017	47.20.017	47.20.117	47.20.317	48.20.117	
37,5	40.12.375	R L	2	47.10.018	47.20.018	47.20.118	47.20.318	48.20.118	
38	40.12.380	R L	2	47.10.019	47.20.019	47.20.119	47.20.319	48.20.119	
38,5	40.12.385	R L	2	47.10.020	47.20.020	47.20.120	47.20.320	48.20.120	
39	40.12.390	R L	2	47.10.021	47.20.021	-	-	-	
39,5	40.12.395	R L	2	47.10.022	47.20.022	-	-	-	
40	40.12.400	R L	2	47.10.023	47.20.023	-	-	-	
40,5	40.12.405	R L	2	47.10.024	47.20.024	-	-	-	
41	40.12.410	R L	2	47.10.025	47.20.025	-	-	-	
41,5	40.12.415	R L	2	47.10.026	47.20.026	-	-	-	
42	40.12.420	R L	2	47.10.027	47.20.027	-	-	-	

ANATOMICAL BAND with tube Upper				TUBES					
				SIMPLE	DOUBLE		TRIPLE	TRIPLE+ LINGUAL	
Prescription				ROTH	ROTH	M.B.T. *L.B.	ROTH	ROTH	LINGUAL
Torque				-14°	-14°	-14°	-14°	-14°	-
Offset				14°	14°	10°	14°	14°	-
Tube ø				-	-	.045"	.045"	.045"	-
Slot				.022"	.022"	.022"	.022"	.022"	2x .036"
									
Perimeter mm	Bands	R / L	Qty	20.10.211 / 20.10.212	20.10.417 / 20.10.418	20.35.402 / 20.35.403	20.10.611 / 20.10.612	20.10.611 / 20.10.612	20.61.206 / -
32	40.22.320	R L	2	47.11.007	-	-	-	-	
32,5	40.22.325	R L	2	47.11.008	-	-	-	-	
33	40.22.330	R L	2	47.11.009	-	-	47.31.009	48.31.109	
33,5	40.22.335	R L	2	47.11.010	-	-	47.31.010	48.31.110	
34	40.22.340	R L	2	47.11.011	47.21.011	47.21.111	47.31.011	48.31.111	
34,5	40.22.345	R L	2	47.11.012	47.21.012	47.21.112	47.31.012	48.31.112	
35	40.22.350	R L	2	47.11.013	47.21.013	47.21.113	47.31.013	48.31.113	
35,5	40.22.355	R L	2	47.11.014	47.21.014	47.21.114	47.31.014	48.31.114	
36	40.22.360	R L	2	47.11.015	47.21.015	47.21.115	47.31.015	48.31.115	
36,5	40.22.365	R L	2	47.11.016	47.21.016	47.21.116	47.31.016	48.31.116	
37	40.22.370	R L	2	47.11.017	47.21.017	47.21.117	47.31.017	48.31.117	
37,5	40.22.375	R L	2	47.11.018	47.21.018	47.21.118	47.31.018	48.31.118	
38	40.22.380	R L	2	47.11.019	47.21.019	47.21.119	47.31.019	48.31.119	
38,5	40.22.385	R L	2	47.11.020	47.21.020	47.21.120	47.31.020	48.31.120	
39	40.22.390	R L	2	47.11.021	47.21.021	47.21.121	47.31.021	48.31.121	
39,5	40.22.395	R L	2	47.11.022	-	-	47.31.022	48.31.122	
40	40.22.400	R L	2	47.11.023	-	-	47.31.023	48.31.123	
40,5	40.22.405	R L	2	47.11.024	-	-	47.31.024	48.31.124	
41	40.22.410	R L	2	47.11.025	-	-	47.31.025	48.31.125	
41,5	40.22.415	R L	2	47.11.026	-	-	47.31.026	48.31.126	
42	40.22.420	R L	2	47.11.027	-	-	47.31.027	48.31.127	

NOTE:
*L.B. = Tube for LIP BUMPER
(Active Lip Plate)

Universal band kits

Size	Qty	LOWER REF 7 6 | 6 7	Upper REF 7 6 | 6 7	Color
	260	-	40.02.900	
	260	40.01.900	-	
	130	-	40.02.901	
	130	40.01.901	-	

Kits with Bands

In practical packaging, the band kits contain all the band sizes, in quantities sized according to their use. Standardized in red for the upper bands packages and blue for the lower ones. The boxes are also sold individually. (See page .124).



Kits composition

Perimeter mm	28	28,5	29	29,5	30	30,5	31	31,5	32	32,5	33	33,5	34	34,5	35	35,5	36	36,5	37	37,5	38	38,5	39	39,5	40	40,5	41	41,5	42
Qtd / Qty	4	4	4	4	4	4	4	8	8	8	8	8	10	10	10	14	16	16	20	20	20	18	10	8	4	4	4	4	4

Pack of 260 Units.

Kits composition

Perimeter mm	28	28,5	29	29,5	30	30,5	31	31,5	32	32,5	33	33,5	34	34,5	35	35,5	36	36,5	37	37,5	38	38,5	39	39,5	40	40,5	41	41,5	42
Qtd / Qty	2	2	2	2	2	2	2	4	4	4	4	4	5	5	5	7	8	8	10	10	9	5	4	2	2	2	2	2	2

Pack of 130 Units.




Blade thickness x height Lâmina x Altura	Qty	REF
0,10 x 3,2mm .004" x .125"	8 metros	40.91.001
0,10 x 3,8mm .004" x .150"	8 metros	40.91.002
0,10 x 4,5mm .004" x .180"	8 metros	40.91.003
0,15 x 4,5mm .006" x .180"	8 metros	40.91.004

Rolls of Banding Tape

The banding tapes are made of polished, malleable stainless steel and are used for bandaging canines and other teeth.


Archs




Full-size diagram of the Morelli intraoral archwires

CrNi Diagrams

Cooper and Stainless Steel



Lower




Upper

NiTi Diagrams


Superelastic, Thermo-plus, Rhodium and Reverse Curve

We offer three sizes of archwires, for better adaptation to different arch sizes and clinical objectives.




Small

For smaller arches, in general lower



Medium

For middle, lower or upper arches



Large

For larger arches, generally higher

Intraoral Archwire CrNi

Size	Qty	Upper	Lower
Identification			
.012”	10	50.60.000	50.70.000
.014”	10	50.60.001	50.70.001
.016”	10	50.60.002	50.70.002
.018”	10	50.60.003	50.70.003
.020”	10	50.60.004	50.70.004
.016” x .016”	10	50.61.001	50.71.001
.016” x .022”	10	50.62.001	50.72.001
.017” x .025”	10	50.62.002	50.72.002
.018” x .025”	10	50.62.003	50.72.003
.019” x .025”	10	50.62.004	50.72.004
.021” x .025”	10	50.62.005	50.72.005




CrNi Archwire

Morelli stainless steel orthodontic archwires have three types of cross-section: round, square and rectangular.



Intraoral Archwire Superelastic NiTi

Medidas *Size*	Qtd *Qty*	Grande *Large*	Médio *Medium*	Pequeno *Small*
Identificação / *Identification*				
.010”	10	-	50.70.010	-
.012”	10	50.60.011	50.70.011	50.80.011
.014”	10	50.60.012	50.70.012	50.80.012
.016”	10	50.60.013	50.70.013	50.80.013
.018”	10	50.60.014	50.70.014	50.80.014
.020”	10	50.60.015	50.70.015	50.80.015
.016” x .016”	10	50.61.002	50.71.002	50.81.002
.016” x .022”	10	50.62.011	50.72.011	50.82.011
.017” x .025”	10	50.62.012	50.72.012	50.82.012
.018” x .025”	10	50.62.013	50.72.013	50.82.013
.019” x .025”	10	50.62.014	50.72.014	50.82.014
.021” x .025”	10	50.62.015	50.72.015	50.82.015



NiTi Superelastic Archwire

NiTi (nickel-titanium) alloys are relevant in orthodontic biomechanics because they have superelasticity and shape memory effects. They release light and constant forces compared to stainless steel. They are not permanently deformed, i.e. the wires cannot be bent. Morelli NiTi archwires come in Small, Medium and Large models.



Superelastic Archwire NiTi Reverse Curve

These archwires are useful in situations where the presence of a curve in the arch is important for orthodontic biomechanics. They are available in different thicknesses and cross-sections (round, square and rectangular).



Intraoral Archwire Superelastic NiTi Reverse Curve - Spee

Size		Qty	Large	Medium
Identification				
◯	.012”	10	50.60.021	50.70.021
◯	.014”	10	50.60.022	50.70.022
◯	.016”	10	50.60.023	50.70.023
◯	.018”	10	50.60.024	50.70.024
◯	.020”	10	50.60.034	50.70.034
□	.016” x .016”	10	50.61.003	50.71.003
▭	.016” x .022”	10	50.62.021	50.72.021
▭	.017” x .025”	10	50.62.022	50.72.022
▭	.019” x .025”	10	50.62.023	50.72.023



Archwire Superelastic NiTi


Made from Nickel-TiTanium (Nitinol) and covered with high-tech polymer. Morelli aesthetic intraoral archwires have a color close to tooth enamel, a resistant covering, flexible and biocompatible material, resistance to staining and color change.



Intraoral Archwire Superelastic NiTi Aesthetic


Size		Qty	Large	Medium
Identification				
◯	.012”	05	50.60.035	50.70.035
◯	.014”	05	50.60.036	50.70.036
◯	.016”	05	50.60.037	50.70.037
◯	.018”	05	50.60.038	50.70.038
▭	.017” x .025”	05	50.62.036	50.72.036
▭	.019” x .025”	05	50.62.037	50.72.037

*Aesthetic arcs equal to and above 0.16 apply only to slot 22



Intraoral Archwire NiTi Thermo Plus

Size		Qty	Large	Medium	Small
Identification					
◯	.012”	10	50.60.224	50.70.224	50.80.224
◯	.014”	10	50.60.225	50.70.225	50.80.225
◯	.016”	10	50.60.226	50.70.226	50.80.226
◯	.018”	10	50.60.227	50.70.227	50.80.227
◯	.020”	10	50.60.228	50.70.228	50.80.228
□	.016” x .016”	10	50.61.204	50.71.204	50.81.204
▭	.016” x .022”	10	50.62.224	50.72.224	50.82.224
▭	.017” x .025”	10	50.62.225	50.72.225	50.82.225
▭	.018” x .025”	10	50.62.226	50.72.226	50.82.226
▭	.019” x .025”	10	50.62.227	50.72.227	50.82.227
▭	.021” x .025”	10	50.62.228	50.72.228	50.82.228




Thermo-Plus Archwire

Thermo-Plus heat-activated NiTi alloy undergoes microstructural transformation with thermal change. When cooled, the archwire becomes malleable and inactive. When it reaches the temperature of the oral environment (above 27°c), it activates and acquires shape memory, generating orthodontic force for tooth movements.



Intraoral Archwire Cooper NiTi

Size		Qty	Upper	Lower
Identification				
◯	.014”	10	50.60.150	50.70.150
◯	.016”	10	50.60.151	50.70.151
◯	.018”	10	50.60.152	50.70.152
▭	.014” x .025”	10	50.62.150	50.72.150
▭	.016” x .022”	10	50.62.151	50.72.151
▭	.017” x .025”	10	50.62.152	50.72.152
▭	.018” x .025”	10	50.62.153	50.72.153
▭	.019” x .025”	10	50.62.154	50.72.154



Copper archwire NiTi

Morelli Copper NiTi intraoral archwires have been specially designed to provide the best clinical performance for the self-ligating appliance system. The addition of 6% copper to the alloy's composition enhances the nickel-titanium thermal properties, resulting in the delivery of optimum forces, with unrivaled efficiency of tooth movement.



Intraoral Archwire Rhodium NiTi

Developed from Superelastic NiTi archwires, rhodium-plated archwires provide excellent aesthetics combined with mechanical performance, maintaining the release of light and constant clinical forces throughout treatment. Morelli Rhodium Arches are coated with a thin layer of rhodium, a precious metal with a silver-white color, which gives the metal a whitish appearance, resembling the color of teeth and guaranteeing an excellent visual result.



Intraoral Archwire Cooper Superelastic NiTi Rhodium

Size		Qty	Large	Medium
Identification				
◯	.012”	01	50.60.135	50.70.135
◯	.014”	01	50.60.136	50.70.136
◯	.016”	01	50.60.137	50.70.137
◯	.018”	01	50.60.138	50.70.138
▭	.017” x .025”	01	50.62.136	50.72.136
▭	.019” x .025”	01	50.62.137	50.72.137



Intraoral Archwires Kit for self-ligating brackets

The packaging itself is a booklet explaining the archwires application and their respective forces. The archwires and suggestions for use presented here do not dispense with the need for careful assessment by a qualified professional, in order to adapt the treatment to the diagnosis and biological response of each patient.



Intraoral Archwire kit for Self-ligating

Type	Qty	Size	REF
Superelastic NiTi	20 un*	Ø.010”	50.14.901
Thermo-Plus NiTi		Ø.012” / Ø.014” / Ø.016” .016”x .016” / .016”x .022” .017”x .025” / .021”x .025”	
Chromo Níquel CrNi		.019” x .025”	
Titanium - Molybdenum TMA		.017” x .025”	

***10 pairs of different archwires**
Each size contains 2 archwires, one upper and one lower.

****10 different pairs***
Each size contains 2 archwires, lower and upper.

SSEQUÊNCIA DE ARCOS MOST COMPLETE ARCHWIRES KIT ON THE MARKET

Intraoral Archwires Kit for Self-ligating

Type	Qty	Size	REF
Superelastic NiTi	18 un*	Ø.012" / Ø.014" Ø.016" / Ø.018" .016" x .025" .017 x .025" .019 x .025"	50.14.900
Chromo Níquel CrNi		.019 x .025"	
Titanium - Molybdenum TMA		.017 x .025"	

***9 pairs of different archwires**
Each size contains 2 archwires, one upper and one lower.

****9 different pairs***
Each size contains 2 archwires, lower and upper.

Intraoral Archwires Kit for Self-ligating

Kit composed of NiTi, CrNi and TMA archwires, with arch sizes to be used with conventional brackets.



Archwire kit for Self-ligating Aesthetic Brackets*A*

Type	Qty	Size	REF
Superelastic NiTi Aesthetic	12 un*	Ø.012" / Ø.014" Ø.016" / Ø.018"	50.14.902
		.017" x .025" .019" x .025"	

***6 pairs of different archwires**
Each size contains 2 archwires, one upper and one lower.

****6 different pairs***
Each size contains 2 archwires, lower and upper.

Archwire kit for Self-ligating Aesthetic Brackets

Sequence composed of aesthetic superelastic archwires NiTi, in a kit of archwire sizes to be used with the Aesthetic Brackets.



Extraoral Arches

Morelli extraoral archwires are made of stainless steel, which is highly resistant and malleable. Together with the various types of cervical caps and splints, these devices offer versatility in the angle of application, allowing orthodontic or orthopedic forces to be obtained. It comes with an exclusive gray bag for archwires.





Extraoral arches

Type	Size	Qty	REF
Standard	 Ø1,15mm (.045")	1	50.10.001
	 Ø1,30mm (.051")	1	50.11.001

Type	Size	Qty	REF
Medium Hard	 Ø1,15mm (.045")	1	50.10.002
	 Ø1,30mm (.051")	1	50.11.002
Medium Hard Small	 Ø1,15mm (.045")	1	50.10.003
Medium White	 Ø1,15mm (.045")	1	50.10.005



It accompanies an exclusive gray bag
Accompanies exclusive grey bag

Lip Bumper

Morelli active lip plates are made of stainless steel and injected plastic, with excellent polish and surface smoothness. They can be used to reinforce anchorage, move molars and treat labial hypotonia.



REF 1: Limp bumper without hook



REF 2: Limp bumper with hook

Lip bumper


Colors	Size	Qty	REF 1	REF 2
Blue	Ø 1,15mm - .045"	1	50.40.002	50.50.002
Pearl	Ø 1,15mm - .045"	1	50.40.006	50.50.006
Pink	Ø 1,15mm - .045"	1	50.40.007	50.50.007

Wires



Wire in Stick Morest CrNi

Stainless steel wires in stick, available in the diameters required for all stages of orthodontic treatment.




Wire in stick Morest

Size	Qty	Code	Qty	REF
○ .012" - 0,30mm	10	55.01.012	-	-
○ .014" - 0,35mm	10	55.01.014	-	-
○ .016" - 0,40mm	10	55.01.016	-	-
○ .018" - 0,45mm	10	55.01.018	-	-
○ .020" - 0,50mm	10	55.01.020	-	-
○ .024" - 0,60mm	10	55.01.024	20	55.06.060
○ .028" - 0,70mm	10	55.01.028	20	55.06.070
○ .032" - 0,80mm	10	55.01.032	20	55.06.080
○ .036" - 0,90mm	10	55.01.036	20	55.06.090
○ .040" - 1,00mm	10	55.06.100	-	-
○ .047" - 1,20mm	10	55.06.120	-	-
□ .016" x .016" 0,40 x 0,40mm	10	55.02.002	-	-
□ .016" x .022" 0,40 x 0,55mm	10	55.03.011	-	-
□ .017" x .022" 0,43 x 0,55mm	10	55.03.016	-	-
□ .017" x .025" 0,43 x 0,63mm	10	55.03.012	-	-
□ .018" x .025" 0,45 x 0,63mm	10	55.03.013	-	-
□ .019" x .025" 0,48 x 0,63mm	10	55.03.014	-	-
□ .021" x .025" 0,53 x 0,63mm	10	55.03.015	-	-

TwistFlex Wire in Stick CrNi

Twist-Flex wire is indicated for the first stages of orthodontic treatment, as well as for making lingual retainer bars.



Twist-flex CrNi Wire in stick

	Size	Qty	REF
3 fios	.015" - 0,38mm	10	55.04.315
3 fios	.018" - 0,45mm	10	55.04.318
3 fios	.020" - 0,50mm	10	55.04.320
4 fios	.015" - 0,38mm	10	55.04.415



Twist-flex CrNi wire in stick for lingual bar

	Size	Qty	REF
3 fios	.032" - 0,80mm	10	55.04.332


#DICA Also check out the TwistForm wire for retention on page: 53

Wire in Stick Blue tip morealloy CoCr

Cobalt-chromium wires, with strength levels suitable for the bioprogressive technique

Blue Tip morealloy CoCr Wire in stick

Size	Qty	REF
.016" x .016" 0,40 x 0,40mm	10	55.02.001
.016" x .022" 0,40 x 0,55mm	10	55.03.001
.017" x .025" 0,43 x 0,63mm	10	55.03.002




TMA Wire in Stick

The special titanium-molybdenum alloy of Morelli TMA wires brings numerous advantages to the completion of orthodontic treatments due to their bending, welding and shape memory characteristics, as well as the high elasticity of the wires.

The TMA brand is the exclusive property of Dental Morelli Ltda. in Brazil

TMA Wire in stick

Size	Qty	REF
.012" - 0,30mm	10	56.01.001
.014" - 0,35mm	10	56.01.002
.016" - 0,40mm	10	56.01.003
.018" - 0,45mm	10	56.01.004
.020" - 0,50mm	10	56.01.005
.032" - 0,80mm	10	56.01.006
.036" - 0,90mm	10	56.01.007
.016" x .016" 0,40 x 0,40mm	10	56.02.001
.016" x .016" 0,40 x 0,55mm	10	56.03.001
.017" x .025" 0,43 x 0,63mm	10	56.03.003
.018" x .025" 0,45 x 0,63mm	10	56.03.004
.019" x .025" 0,48 x 0,63mm	10	56.03.005
.021" x .025" 0,53 x 0,63mm	10	56.03.006



Wires in Rolls

Stainless steel wires, with a round cross-section, are used for various purposes in orthodontics, such as making removable appliances. Available in 50 g or 500 g rolls.

Hard elastic wire CrNi - 50 g

Medidas *Size*	Metragem Aprox. *Aproximated measure*	REF
.012" - 0,30mm	90m	55.01.030
.014" - 0,35mm	66m	55.01.035
.016" - 0,40mm	50m	55.01.040
.018" - 0,45mm	40m	55.01.045
.020" - 0,50mm	32m	55.01.050
.024" - 0,60mm	22m	55.01.060
.028" - 0,70mm	16m	55.01.070
.032" - 0,80mm	12m	55.01.080
.036" - 0,90mm	10m	55.01.090
.040" - 1,00mm	08m	55.01.100
.047" - 1,20mm	05m	55.01.120



CrNi Wire in stick

Stainless steel wires, available in the diameters required for all stages of orthodontic treatment.



Hard elastic wire CrNi - 500 g

Size	Approximated measure	REF
.028" x .070mm	165m	55.01.570
.032" x .080mm	126m	55.01.580
.036" x .090mm	100m	55.01.590
.040" x 1,00mm	81m	55.01.591



Hard elastic half-cane wire CrNi

Size	Approximated measure	REF
Ø 2,00 x 1,00mm	5m	55.01.200
Ø 1,30 x 0,65mm	5m	55.01.130
Ø 1,50 x 0,75mm	5m	55.01.150
Ø 1,75 x 0,90mm	5m	55.01.175




Wire for ligature CrNi - 50g

Size	Approximated measure	REF
.008" - 0,20mm	201m	55.01.208
.010" - 0,25mm	130m	55.01.210
.012" - 0,30mm	90m	55.01.212




Aesthetic wire for CrNi ligature in roll

Size	Approximated measure	REF
(D1) Internal Wire .010" - 0,25mm (D2) External Wire .014" - 0,35mm	10m	55.01.220





Elastics





Elastic ligatures
with antimicrobial additive

The new line of ligatures features new, super modern and exclusive colors for **modulars and sticks**, all with an antimicrobial additive, an innovation to make the elastics more durable and durable.

Antimicrobial Technology:


The incorporation of the antimicrobial additive into the polymer matrix prevents the development and growth of fungi and bacteria. This innovation **reduction of stains and unwanted odors, as well as preserving the elastics mechanical properties**.

Exclusive Colors:

The ligatures have carefully selected, top-quality pigments to create a more stable, modern and vibrant color palette. **There are 9 exclusive options** to meet different tastes and styles.

Shape memory and mechanical properties:

By using **excellent raw materials** and the antimicrobial additive, we produce elastics with better mechanical properties. In addition, we heat-treat the material, ensuring greater shape memory for the elastics.



Plic Tweezers
for modular

Perfectly fitting tip, for easily applying elastic in the bracket.

Also great for applying stops, arches, springs and hooks.

REF 75.01.028

Elastic ligatures

Colors	Stick		Modular	
Amethyst Purple	1000 un.	60.03.166	500 un.	60.06.166
Bubblegum Pink	1000 un.	60.03.168	500 un.	60.06.168
Dragon Fruit Pink	1000 un.	60.03.165	500 un.	60.06.165
Coral Orange	1000 un.	60.03.167	500 un.	60.06.167
Lime Green	1000 un.	60.03.161	500 un.	60.06.161
Caribbean Blue	1000 un.	60.03.162	500 un.	60.06.162
Sky Blue	1000 un.	60.03.163	500 un.	60.06.163
Teal Blue	1000 un.	60.03.164	500 un.	60.06.164
Graphite Gray	1000 un.	60.03.160	500 un.	60.06.160

Printed colors are approximate and may change.

Request your color display here

Elastic ligatures

Colors	Qty. Rngs	Stick	Modular
Cinza	5000	60.03.701	-
Cristal	5000	60.03.700	-
Amarelo	1000	60.03.305	60.06.104
Azul	1000	60.03.308	60.06.107
Azul Bebê	1000	60.03.304	60.06.103
Azul Cristal	1000	60.03.354	60.06.154
Azul Marinho	1000	60.03.316	60.06.115
Branco Renda	1000	60.03.327	60.06.127
Cristal	1000	60.03.300	60.06.100
Cinza	1000	60.03.301	60.06.101
Laranja Cristal	1000	60.03.352	60.06.152
Lilás	1000	60.03.328	60.06.128
Marfim	1000	60.03.320	60.06.120
Pink Cristal	1000	60.03.355	60.06.155
Prateado	1000	60.03.317	60.06.117
Preto	1000	60.03.311	60.06.110
Rosa	1000	60.03.312	60.06.111
Rosa Bebê	1000	60.03.319	60.06.119
Rosa Pérola	1000	60.03.326	60.06.126
Roxo	1000	60.03.313	60.06.112
Sortida	1000	60.03.302	60.06.118
Sortida Cristal	1000	60.03.358	60.06.158
Uva Cristal	1000	60.03.359	60.06.159
Verde Cristal	1000	60.03.353	60.06.153
Verde limão Cristal	1000	60.03.350	60.06.150
Verde Musgo	1000	60.03.307	60.06.106
Verde Mar	1000	60.03.315	60.06.114
Verde Pérola	1000	60.03.324	60.06.124
Vermelho	1000	60.03.310	60.06.109
Vermelho Cristal	1000	60.03.351	60.06.151
Vinho	1000	60.03.322	60.06.122

Printed colors are approximate and may change


Elastic for ligature

Morelli elastics for ligature use polyurethane controlled in their own manufacturing process and non-toxic pigmentation. It comes in different colors, in stick and modular versions. These elastics maintain efficient memory levels and low deformation rates in the oral environment.

Modular Stick




Modular





Request your color display here

Elastic chain

Elastic chain is useful for closing spaces, retracting anterior teeth and pulling teeth. They are made from high-quality material, offering elasticity and memory. Available in short, medium and long link versions, colored and in packs of 4.5 m (gray and crystal) and 1.5 m (other colors).







Elastic chain

Color	Lenght	Short	Medium	Long
Grey	4,5m	60.05.501	60.05.511	60.05.521
Crystal	4,5m	60.05.500	60.05.510	60.05.520
Yellow	1,5m	60.05.105	60.05.205	60.05.305
Blue	1,5m	60.05.108	60.05.208	60.05.308
Baby BlueBebê	1,5m	60.05.104	60.05.204	60.05.304
Crystal blue	1,5m	60.05.154	60.05.254	60.05.354
Navy Blue	1,5m	60.05.116	60.05.216	60.05.316
Lace White	1,5m	60.05.127	60.05.227	60.05.327
Crystal Orange	1,5m	60.05.152	60.05.252	60.05.352
Lilac	1,5m	60.05.128	60.05.228	60.05.328
Ivory	1,5m	60.05.120	60.05.220	60.05.320
Chrystal Pink	1,5m	60.05.155	60.05.255	60.05.355
Silver	1,5m	60.05.117	60.05.217	60.05.317
Black	1,5m	60.05.111	60.05.211	60.05.311
Rose	1,5m	60.05.112	60.05.212	60.05.312
Baby Pink	1,5m	60.05.119	60.05.219	60.05.319
Pearl Rose	1,5m	60.05.126	60.05.226	60.05.326
Purple	1,5m	60.05.113	60.05.213	60.05.313
Crystal Grape	1,5m	60.05.159	60.05.259	60.05.359
Crystal Green	1,5m	60.05.153	60.05.253	60.05.353
Crystal Lime Green	1,5m	60.05.150	60.05.250	60.05.350
Moss Green	1,5m	60.05.107	60.05.207	60.05.307
Sea Green	1,5m	60.05.115	60.05.215	60.05.315
Pearl Green	1,5m	60.05.124	60.05.224	60.05.324
Red	1,5m	60.05.110	60.05.210	60.05.310
Crystal Red	1,5m	60.05.151	60.05.251	60.05.351
Wine	1,5m	60.05.122	60.05.222	60.05.322

Printed colors are approximate and may change



Intraoral Elastics

Black LATEX

The Morelli Intraoral Elastics line now has a Black Latex version with an antimicrobial additive.

An innovation to bring more aesthetics and durability to elastics.


Black LATEX


LATEX


SILICONE

More discretion

The lighter tone of natural latex and silicone stands out in the patient's oral cavity. By using black latex, we have more discretion when opening the mouth, which can help with acceptance during treatment.

Antimicrobial Technology

The incorporation of the antimicrobial additive into the elastomer matrix prevents the development and attachment of fungi and bacteria. This innovation preserves the elastics mechanical properties and ensures that the color of the Black Latex is maintained.

Size	Code	Quantity
Medium 1/8" = 3,20mm	60.02.410	1000
Medium 3/16" = 4,80mm	60.02.411	1000
Medium 1/4" = 6,35mm	60.02.412	1000
Medium 5/16" = 7,94mm	60.02.413	1000



Elastics ligatures

Color	Qty.Rings	Modular
Pearl Blue	500	60.06.010

Pearl Blue Elastic

Pearl Blue Morelli elastics are ideal for all orthodontic treatments, and especially suitable for treatment with aesthetic brackets. Due to its high resistance to staining and bluish color, it has an excellent aesthetic result!



Elastics chain

Color	Lenght	Short	Medium	Long
Pearl Blue	1,5m	60.05.101	60.05.201	60.05.301

Latex intraoral elastic

The natural latex elastic properties are widely used in orthodontics for dental intercuspations of various levels. They have a min scent.



3/4" Latex intraoral elastics

Force	Color	Size	Qty	REF
Light	Natural	3/4” = 19,05mm	100	60.01.330
Medium		3/4” = 19,05mm	100	60.01.340

Protocols and clinical tips on the use of intraoral elastics


Latex intraoral elastics

Force	Colors	Size	Qty	REF
Light	Natural	1/8” = 3,20mm	1000	60.01.200
		3/16” = 4,80mm	1000	60.01.201
		1/4” = 6,35mm	1000	60.01.302
		5/16” = 7,94mm	1000	60.01.303
		3/8” = 9,53mm	1000	60.01.304
Medium	Natural	1/8” = 3,20mm	1000	60.01.310
		3/16” = 4,80mm	1000	60.01.311
		1/4” = 6,35mm	1000	60.01.203
		5/16” = 7,94mm	1000	60.01.205
		3/8” = 9,53mm	1000	60.01.314
Heavy	Natural	1/8” = 3,20mm	1000	60.01.320
		3/16” = 4,80mm	1000	60.01.321
		1/4” = 6,35mm	1000	60.01.322
		5/16” = 7,94mm	1000	60.01.323
		3/8” = 9,53mm	1000	60.01.208



Latex extraoral elastics

Force	Color	Size	Qty	REF
Light	Natural	1/2" = 12,7mm	100	60.01.210
Medium		1/2" = 12,7mm	100	60.01.211
Heavy		1/2" = 12,7mm	100	60.01.212



Silicone Intraoral Elastic

Morelli Silicone intraoral elastics have been developed to meet the mechanical needs demanded by brackets and accessories with an aesthetic, biocompatible and anti-allergenic appeal. They have a min scent.

Silicone intraoral elastics

Force	Color	Size	Qty	REF
Medium	Translucent	1/8" = 3,20mm	1000	60.08.110
		3/16" = 4,80mm	1000	60.08.111
		1/4" = 6,35mm	1000	60.08.112
		5/16" = 7,94mm	1000	60.08.113

The nominal force will stabilize only after the first elongation.



Separator elastic


Morelli separator elastics meet the specifications for elasticity and resistance, promoting an efficient opening of the interproximal space.

Elastic separator

Force	Color	Size	Qty	REF
In Bulk	Blue	5/32" = 4,0mm	1000	60.04.200
	Green	3/16" = 4,80mm	1000	60.04.400



Modular	Blue	5/32" = 4,0mm	500	60.04.201
	Green	3/16" = 4,80mm	500	60.04.401



Protection Tube

This tube has the function of protecting the mucosa and the wires themselves, providing reinforcement and comfort for the patient.

In addition to protection, the colored tubes provide a decorative effect and the possibility of combining colors with the other ligatures used in orthodontic treatment, promoting the aesthetic appeal that the patient desires in a safe way.



Tube for protection

Color	Size	Diameter	REF
Crystal	3,0m	0.75mm	60.05.400
Grey	3,0m	0.75mm	60.05.401
Blue	3,0m	0.75mm	60.05.403
Crystal Pink	3,0m	0.75mm	60.05.405
Crystal	3,0m	0.95mm	60.05.410
Grey	3,0m	0.95mm	60.05.411
Blue	3,0m	0.95mm	60.05.402
Crystal Pink	3,0m	0.95mm	60.05.404



***Diameter: 0.75 mm for .012" to .016"x.025" wires 0.95 mm for .017"x.02 5" to .021"x.02 5" wires.**


Rotator Elastic

Rotator elastics are positioned on the bracket fins to promote rotational movement of the tooth



Rotator elastics


Color	Qty	REF
Crystal	100	60.07.100
Grey	100	60.07.101






Universal Rotator Elastics

Color	Qty	REF
Crystal	100	60.07.200
Grey	100	60.07.201



Check out how to use


Functional Orthopedics



ortho protector

MULTIPURPOSE
RAISING PLATE

The bite raising plate with no moldings and instant adaptation to any patient



Ortho Protector is a bite raising plate that acts as a muscle relaxant.

Temporarily treats occlusal stress problems such as **bruxism** and temporomandibular disorders. Practical and comfortable, even for patients with braces.

It also helps protect laminate veneers and rehabilitate after implant surgery.

Resistant material

Made of polyurethane, which provides durability, even with occlusal tightening.

Occlusal plates

They easily fit all types and sizes of arches, including children from 5 years old.

Front strap

It fits anatomically on the buccal side of the lower incisors.



Size adjustments

There are 4 positions available for better adaptation and comfort during use.

What are the indications for Ortho Protector?



Bruxism

The Ortho Protector can be used to control bruxism, even in children over 5 years of age, as it is fully adjustable to the arches.

Occlusal Stress

Ortho Protector provides occlusal isolation, which relieves the neural system, inhibiting the development of muscle tension and orofacial pain.

Temporomandibular Disorders

Indicated as a resource for immediate relief of orofacial pain caused by Temporomandibular Disorders (TMD)

Orthodontics

The Ortho Protector can be used by patients during orthodontic treatment, as its design is fully adjustable and increases orthodontic movement by providing occlusal raising. Also indicated in cases of maxillary disjunction and traction, as it helps to dislodge the bite. It also provides muscle relaxation and relief from the discomfort caused after arch and elastic maintenance.

Laminate Veneers

After the installation of laminate veneers or “dental contact lenses”, Ortho Protector is a safe option for nighttime protection, reducing the risk of cracks or fractures.

Prosthetic and Implant Rehabilitation

The use of Ortho Protector is also highly recommended for the nighttime protection of patients undergoing rehabilitation after the installation of removable or implanted prostheses.

Harmful Habits

The use of Ortho Protector can be indicated for patients who have habits that damage parts of the oral cavity, such as the habit of biting the buccal mucosa (inner part of the cheeks).

Watch the video with a full explanation of how to use Ortho Protector.

Color	Qty.	REF
Blue	1	65.03.001
Pink	1	65.03.002

It comes with a box for storage and transportation, and instructions for use for the patient.



Mini Hyrax Expansion screw

A new treatment option for mandibular atresia

The Mini Hyrax Expansion screw, also known as the “Williams Expansion screw”, is an effective alternative for treating severe crowding of the lower teeth by expanding the lower arch.

Unlike the Hyrax Expansion screw, which allows orthopedic movement by opening the palate, the Mini Hyrax Expansion screw only promotes dentoalveolar movement, for transverse expansion of the lower arch.

Openings


Available in two opening sizes: 9 mm and 11 mm. Mini Hyrax Expansion screws are activated by means of a central screw.

Stainless steel

Made entirely of stainless steel, guaranteeing high resistance.

Small size

It can be installed on the lingual side of the anterior teeth, with the stems anchored in bands, for tooth movements in the mandible.



Anti-return

Brake bushing that prevents the screw from returning unintentionally, maintaining activation.







Applications

The Mini Hyrax Expansion screw is an effective alternative for expanding the lower arch. It can also be applied to the maxilla to perform palatal expansion, if the orthodontist so chooses.




The Mini Hyrax Expansion screws, like the Palatal Disjunction Expansion screws, are also supplied with a **Long Stem Wrench**.





Opening	Qty.	REF
9 mm	1	65.05.111
11 mm	1	65.05.112

Hyrax Expansion screw Hyrax




Hyrax expansion screws are indicated for correcting deficient maxillomandibular relationships.

The most common clinical application situations are:

- Cases of early mixed/permanent dentition;
- Posterior-superior bi-lateral crossbite as a result of deficiency of the maxilla apical base;
- An anterior or posterior crossbite that is not too severe.

Hyrax - Palatal expansion screw




Opening	Qty	REF
7mm	1	65.05.010
9mm	1	65.05.011
11mm	1	65.05.012
13mm	1	65.05.013

Fan Type

The rapid fan expansion screw is designed for the upper arch anterior expansion and is mounted in the anterior palatine region. It is commonly used to treat patients with craniofacial anomalies (cleft palate)



Rapid fan expansion screw



Opening	Qty	REF
10mm	1	65.05.014



All expansion screws come with an activation key - code 75.01.038 - page 115

**Expansion: | Opening:



Bimler Shield


Size	Qty	REF
	10	65.01.010

Bimler Shield de Bimler

Device that assembles and joins the lower arches of the Bimler elastic shaper. The Bimler Shield is used in cases of malocclusion with or without crowding, and lack of space for canines.




Equiplan flats


Size	Qty	REF
	10	65.01.040



Standard Expansion Screw



Opening	Qty	REF
7mm	10	65.05.001
7mm	1	65.05.101



11mm	10	65.05.009
11mm	1	65.05.109

Hyrax Expansion screw


Morelli expansion screws are used in fixed or mobile appliances to expand or disjunction dental arches. They are available in various sizes, models and openings.



All Expansion Screws are accompanied by an activation key - code 75.02.026 - page 115


*Expansão: *Expansion* | Abertura: *opening*

Universal Expansion screw




Opening	Qty	REF
7mm	10	65.05.002
7mm	1	65.05.102





Opening	Qty	REF
9mm	10	65.05.003
9mm	1	65.05.103






Opening	Qty	REF
11mm	10	65.05.008
11mm	1	65.05.108



Mini Expansion screw




Opening	Qty	REF
6,5mm	10	65.05.006
6,5mm	1	65.05.106



*Expansion: | Opening:

Mini Expansion screw






Opening	Qty	REF
9mm	10	65.05.007
9mm	1	65.05.107

Mini unilateral expansion screw






Opening	Qty	REF
5mm	1	65.05.156

Unilateral expansion screw






Opening	Qty	REF
6,5mm	1	65.05.151

Fan type expansion screw






Opening	Qty	REF
8mm	1	65.05.160

*Expansion: | Opening:

Headcaps and Tractioners



Headcap

They are used for extraoral anchorage, allowing high tensile forces to be applied.



Headgear for Safety Breakaway Device

Size	Qty	REF
Small	1	70.31.001
Medium	1	70.31.002
Large	1	70.31.003



Headgear - Adjustable traction Interlandi

Size	Qty	REF
One size	1	70.35.001

Chin cup

To contain mandibular growth by applying forces in a parietal direction. It comes complete or just the chin cup.



Chin cup complete - high traction

Size	Qty	REF
One size	1	70.36.001



Chin cup - chin traction device


Size	Qty	REF
Onze size	1	70.36.002

Cervical flap

To apply extra-mouth force with cervical traction, using tractioners.



Universal cervical flap for low traction

Size	Qty	REF
	1	70.20.032

Petit Facial Mask

Anchorage device used for the maxilla reverse traction with posteroanterior advancement. The Morelli petit face mask provides a frontal downward and backward tilting action of the jaw at the same time.


The forehead support (1) and chin rest (2) have foam protection, which increases comfort and encourages the patient to use.

Petit Facial Mask

Size	Qty	REF
Blue	1	70.50.001
Rose	1	70.50.002
Gray	1	70.50.003



It is possible to use Morelli's Petit facial mask in two ways: Fixed or Articulated chin cup.



Fixed Chin CupFixa

Articulated Chin Cup



Breakway with traction device


Qty	Light force *350gf*	Medium force *500gf*	Heavy force *750gf*
2	70.42.102	70.42.202	70.42.302

Tractioner

Unlocks the clip spontaneously if the force applied exceeds the limit.



Check how to use



Parietal Elastic Double Clasp - low traction


Size	Qty	REF
One size	1	70.20.011




Instruments



121°C



Elastic ligature applicator / extractor

Size	Qty	REF
 3,4mm	1	75.01.001
 3,4mm	1	75.01.002

Ligature Applicator

This tool is used for removing and attaching elastic ligatures.


They have special stainless steel double active tips.

Check out how to open the Morelli Self-ligating clip correctly


121°C



Elastic Ligature Applicator

Size	Qty	REF
 3,4mm	1	75.01.003
 3,4mm	1	75.01.004

121°C



Bracket positioner

Size		Qty	REF
 mm	3,5mm	1	75.01.006
	4,0mm		
	4,5mm		
	5,0mm		
	2,5mm	1	75.01.039
	3,0mm		
	5,5mm		
	6,0mm		

Bracket Positioner

Made of stainless steel, the Morelli bracket positioner has an exclusive hinged tip, which makes it easier to position on posterior teeth.

This instrument allows the correct positioning of the bracket slot in relation to the occlusal or incisal edge of the tooth.





incisal do dente.

Band seater

Suitable for seating the band on the tooth.

Stainless steel tips with shapes that favor the adaptation of orthodontic bands. Made of gray nylon, they resist the pressure of bites and the temperatures of autoclave sterilization.

Band seater

Size	Qty	Format	REF
	1	Triangular tip	75.01.011
	1	Inverted triangular tip invertida	75.01.012
	1	Round tip	75.01.013
	1	Square tip	75.01.014

121°C



121°C

Lip Retractor

Morelli lip retractors are useful for optimizing isolation during direct bonding of attachments, as well as when taking intraoral photographs.

Lip retractor

Size	Color	Qty	REF
Small	Amber	2	75.01.033
Large	Amber	2	75.01.043



Size	Color	Qty	REF
Small	Amber	2	75.01.036
Large	Amber	2	75.01.046



Size	Color	Qty	REF
Small	White	1	75.01.133
Large	White	1	75.01.143



Size	Color	Qty	REF
Small	White	1	75.01.136
Large	White	1	75.01.146



Disposable impression trays

Practicality, hygiene and precision for moldings

Morelli's disposable impression trays are manufactured using a fully automated process, from injection to sealing the individualized packaging, without manual contact, guaranteeing maximum hygiene and safety for patients. For professionals, more practicality and savings, eliminating the need for cleaning and sterilization.

Three sizes, identified by color

Efficient retention

The internal TNT (non-woven fabric) coating improves the hold of the molding material, with a more reliable result, reducing distortions and material waste. It is also more comfortable for the patient, as there is less material leakage.

Hygiene and safety

As they are disposable, they do not require cleaning and sterilization. It can be opened in the presence of the patient, which promotes greater care with asepsis for both the patient and the professional.

Translucency

The trays' translucent body allows you to see the distribution of the material and control it during the process.



Disposable impression trays

Embaladas individualmente, mas unidas aos pares (superior e inferior), e prontas para uso imediato.

Disposable impression trays

Description	Color	REF	Qty
Disposable impression tray - Upper/lower pair S	Green	75.01.152	Package of 280 x 380 mm with 10 units (5 pairs)
Disposable impression tray - Upper/lower pair M	Pink	75.01.153	
Disposable impression tray - Upper/lower pair L	Blue	75.01.154	







Conventional sterilizable impression trays

Morelli impression trays have excellent alginate retention due to the design of the grooves at the base, as well as allowing a complete impression of the dental arch due to their greater lateral height. They can be sterilized both in autoclave or chemically.

Impression trays

Material	Color	Size	REF	Qty
Plastic	Gray	01	75.01.051	1 pair
	Yellow	02	75.01.052	
	Pink	03	75.01.053	
	Orange	04	75.01.054	
	Green	05	75.01.055	
	Red	06	75.01.056	
	Blue	07	75.01.057	
	Purple	08	75.01.058	

121ºC



RELATIONSHIP BETWEEN TRAYS - UPPER / LOWER


Conventional sterilizable trays

Size	height		width		Length	
	upper	lower	upper	lower	upper	lower
01	24	17	50	50	44	43
02	26	22	61	62	46	46
03	28	22	66	69	50	46
04	28	23	70	67	53	50
05	30	25	70	69	55	52
06	32	25	73	74	57	56
07	32	26	77	75	61	59
08	32	27	85	78	68	64


Disposable Molds

Tam	height		width		Length	
	upper	lower	upper	lower	upper	lower
P	24	21	71	70	56	53
M	25	22	74	73	59	56
G	26	22	77	76	64	60


upper tray



lower tray



NOTE: All dimensions are external and for reference only and may vary by +/- 2.0mm / equivalents for length and width




Tower to form arches

Size	Qty	REF
.016”	1	75.02.005
.017”		
.018”		
.019”		
.021”		

Tower to form arches

The tower is an instrument used to shape orthodontic arches from square or rectangular wire sticks of various thicknesses.


The Morelli Tower is made of stainless steel with an engineering polymer handle and an excellent finish.



Removable cover tube detacher

Size	Qty	REF
 130mm	1	75.01.024

Detacher blade

Size	Qty	REF
 33mm	5	75.01.025

Detacher

Easily removes the cap from single, double and triple convertible tubes. It can also be used with scalpel blades. Detacher blades are available for replacement.



Tensiometer

Color	Force	Qty	REF
	100gf a 1600gf (cN)	1	75.02.006
	50gf a 500gf (cN)	1	75.02.008
	25gf a 250gf (cN)	1	75.02.009

Tensiometer

Mechanical tensiometer used to measure forces in intraoral and extraoral elastics.

They have two tips for easy positioning.

Orthodontic tweezers Ortodôntica

The orthodontic tweezers have been designed to allow safe and perfect grasping of brackets and accessories.

The end of the handle has an exclusive finish to facilitate the seating and positioning of brackets.

Orthodontic tweezers

121ºC


Size	Qty	REF
	1	75.01.022





Short tweezers for bracket position

121ºC

Size	Qty	REF
	1	75.01.064




Tweezers for tube

121ºC


Size	Qty	REF
	1	75.01.027





Plic Tweezers

121ºC

Size	Qty	REF
	1	75.01.028

Plic Tweezers


The "Plic" tweezers, thanks to its recessed tip, offers proven advantages when inserting rubber bands, springs, round stops, etc.

Suggested application



Bands marker

121ºC

Size		Qty	REF
	3,5mm	1	75.01.007
	4,0mm		
	4,5mm		
	5,0mm		

Band Marker



This instrument is used to mark accessory heights on orthodontic bands.

It has a rotating mechanism with sizes from 3.5 mm to 5.0 mm.



Boone gauge

121ºC

Size		Qty	REF
	3,5mm	1	75.01.023
	4,0mm		
	4,5mm		
	5,0mm		
	2,0mm	1	75.01.029
	2,5mm		
	3,0mm		
	3,5mm		
	4,0mm	1	75.01.030
	4,5mm		
	5,0mm		
	5,0mm		

Boone gauge


Made of stainless steel, this instrument has heat-treated tips. It is used to mark bands and position brackets and accessories.

Bender

The Morelli bender is a practical device to facilitate bending thick wires (Hyrax expander rods, etc.). It has openings and slots for fitting wires, allowing the use of wires up to 01.5 mm of thickness.

Wire and rods bender

121°C

Size	Qty	REF
 130mm	1	75.02.001





Bender Distal

Used to bend and finish the orthodontic wire distally, avoiding discomfort for the patient. It comes in versions with a hole (closed) or a slot (open), as well as an appropriate curvature to make it easier to apply.

Distal arm bending - closed shot

121°C

Size	Qty	REF
.022”	1	75.02.022



Distal arm bending - opened shot

121°C

Size	Qty	REF
.022”	1	75.02.023





Arm bending for wire ligature - curve tip

121ºC

Size	Qty	REF
.022”	1	75.02.074

Arm bending for wire ligature

The arm bending for ligature are useful for bending and accommodating metal and elastic ligatures in the channels of brackets and other orthodontic appliances. They have straight and curved tips for better adjustment and application.



Arm bending for wire ligature - straight / curve tip


121ºC

Size	Qty	REF
.022”	1	75.02.075



Hexagonal key for Gurin nut

121ºC

Size	Qty	REF
	1	75.02.011

Hexagonal key

Hexagon key for fastening the Gurin nut (page 49)

Elastic Positioner

The positioners (applicators) make it easier for the patient to place the elastics as instructed by the professional.

Elastic positioner

Color	Qty	REF
Black	20	75.01.015

See retainer box for elastic #3 on page 123




Template

Suitable for making orthodontic cephalometric tracings. It allows the professional to get a better view of the patient's problems, which allows them to make a diagnosis that is closer to the ideal.

Benvenga cephalometric template (1995)

Size	Qty	REF
190 x 200 mm	1	75.01.020




Expansion screw key


Allows the expansion screw to be activated safely and conveniently.

Activation key for expansion screws

121ºC

Model	Qty	REF
Long stem	5	75.01.038
Short stem	10	75.02.026





Zirconia drill bit

Suitable for safely removing orthodontic resin without damaging or scratching tooth enamel.

Use with a low-rotation contra-angle, removing up to 5 arches on average.

Zirconia low rotation drill bit

121ºC

Size		Qty	REF
Ø 3mm	23 mm	5	75.03.001



EFFICIENCY
The Morelli Drill bit has been specially developed to remove the remaining resin after removing the appliance.

REUSABLE
High resistance to autoclave sterilization cycles and easy to clean grooves, as they are straight.

GREAT FINISHING
The quality of the drill bit causes low vibration, which generates less heat and perfectly smooth surfaces.

DIFFERENTIATED MATERIALA
The only removal drill bit on the market sintered entirely in zirconia

DURABILITYA
The Zirconia Drill bit has excellent cutting durability, resistance and does not damage or contaminate teeth. It is also efficient in either direction of rotation, which prevents errors and increases its service life.

SAFETY
Its multilaminated active tip does not hurt the gums, which allows the resin to be removed in a controlled manner, eliminating the risk of cracking, scratching the enamel or damaging the soft tissue.

Disposable point

Ideal for restoring the shine and smoothness of tooth enamel after using the Morelli Zirconia Drill bit.

Use with a low-rotation contra-angle, removing up to 1 arch.

Disposable point for polishing

Size		Qty	REF
Ø 3mm	23 mm	20	75.03.002



EXCELLENT FINISHING
The polishing point provides an optimum finishing to the enamel.

SAFETY
Its rounded shape is safer for the gum and allows access to areas close to it, eliminating any remaining resin and ensuring an excellent finishing.

MATERIAL QUALITY
Its fiberglass polymer composition is abrasive enough to remove adhesive residue, guaranteeing an effective polish that is safe from damaging dental enamel.

PRACTICALITY
Morelli point can be applied dry and in any direction of rotation. As it is disposable, it does not require sterilization.

Oral Hygiene



A different, comfortable and safe way to take care of your oral health.

Developed with a different design to conventional toothbrushes, **Bioflex® brushes have 393 "rubberized" bristles,** which provide excellent plaque removal, protecting gums and mucous membranes from injury.

High-technology

Its great advantage is its "monobloc" injected bristles which, unlike traditional "tufted" bristles, have fewer indentations, making the brush more hygienic. The soft bristles also preserve restorations, orthodontic appliances, implants and other oral structures.

Escova Dental

BioFlex®



Innovative design

The Bioflex® line has been developed to meet the needs of all patients. With its innovative design and anatomical shape, the brush handle fits well in the hand, making brushing easier.

Healthy and thorough cleaning


On the back of the bristles, the Bioflex® brush also has a tongue scraper, which removes tongue plaque.

Toothbrush Bioflex

Bioflex® brushes bring a different, comfortable and safe way to take care of your oral health. Developed with a different design from conventional toothbrushes, Bioflex® brushes have 393 "rubberized" bristles that provide excellent plaque removal, protecting gums and mucous membranes from injury. Their have soft bristles, which also preserve restorations, orthodontic appliances, implants and other oral structures.

Soft Bioflex® Toothbrush

Color	Qty	REF
White/Pink	01 unit	78.01.200
Blue/Lemon	01 unidade	78.01.201
Pink/Purple	01 unidade	78.01.202
White/Blue	01 unidade	78.01.203



Display

Containing 10 assorted brushes, the display is ideal for showing off at the reception desk.

Bioflex® Toothbrush Soft Assorted on Display

Color	Qty	REF
Assorted	01 kit 10 brushes	78.01.210

Kit composed of:
2 brushes (White/Pink), 3 brushes (Blue/Lemon), 3 brushes (Pink/Purple), 2 brushes (White/Blue).





Interdental Brush

Color	Qty	REF
White	01 kit	78.01.001
Blue	01 kit	78.01.002
Yellow	01 kit	78.01.003
Pink	01 kit	78.01.004

Interdental Brush

Interdental brushes are used to clean between the teeth (interproximal), preventing tartar and periodontal disease. They are very effective for cleaning orthodontic appliances and reach places that are difficult to access with other cleaning tools.
Morelli interdental brushes help with prophylaxis, favoring the progress of orthodontic treatment, preventing stains and cavities on the teeth, as well as gum infections





Interdental brush Refill

Size	Qty	REF
Cylindrical P	1 blister 6 units	78.01.101
Cylindrical M	1 blister 6 units	78.01.102
Cylindrical G	1 blister 6 units	78.01.103

White coated wire: Better bristle fixation and patient comfort
Kit composed of:
3 Refills (S,M,L), 1 case and 1 handle
Kit composrd of:
3 brushes, 1 case and 1 cable



Tongue Cleaner

Size	Qty	REF
Blue and pink	1 pair 2 units	78.02.001
Blue and yellow	1 pair 2 units	78.02.002

Tongue Cleaner

The Morelli tongue cleaner has an anatomical and functional design, suitable for removing coating (responsible for bad breath and plaque formation) with maximum efficiency and care, without damaging the tongue papillae.

Compartments



Boxes

As well as being beautiful, Morelli Pocket-size Boxes are super resistant, practical and hygienic. They have a vent at the top of the lid, which ensures air circulation, keeping both the box and its contents dry at all times.

The Pocket-size Box is suitable for transporting and storing various items, such as: Retention plates, Bruxism plates, Whitening plates, Mobile braces and Aligners.

Pocket-size cases for mobile appliances and elastics

Colors	Qty	REF
Blue/Withe	10	80.04.211
Pink/Withe	10	80.04.212
Blue/Withe	05	80.04.214
Black Piano	05	80.04.215
Turquoise Blue	05	80.04.216
Navy Blue	05	80.04.217
Cherry Red	05	80.04.218
Magenta Rose	05	80.04.219
Tangerine Orange	05	80.04.220
Lime Green	05	80.04.221



Pocket-size cases for mobile appliances and elastics

Colors	Qty	Box	REF
Assorted	10	Removable Appliances	80.03.004



Colors	Qty	Box	REF
Assorted	10	Removable Appliances	80.04.001



Colors	Qty	Box	REF
Assorted	10	For Elastics	80.04.002

Accompany Device for place and remove elastics!

Elastic positioner single code 75.01.015 on page115





Zip bag for elastics

Colors	Qty	Code
78mm / 97mm	50	80.05.001




Organizer Boxes

With a rounded bottom and band and bracket numbering, Morelli boxes are ideal for storing accessories, making your practice more practical and organized.

They are made of high-impact polystyrene.




Box for upper bands

Size	Use	Qty	REF
	Universal	1	80.02.011




Box for lower bands

Size	Use	Qty	REF
	Universal	1	80.02.021




Box for brackets

Size	Qty	REF
	1	80.02.030




Multiuse organizer box

Size	Qty	REF
	1	80.06.020

Accompanied by 64 adhesive labels, which can be sticked in the box to identify the materials.

Factory **Tour**



The Factory Tour Program reflects our philosophy of meeting the needs and suggestions of the market.

That's why our doors are always open to welcome orthodontic students and professionals.

more than 20 years of program

More than 15 thousand visitors

Largest orthodontic center of South America

Come and see Morelli for yourself

SCHEDULE YOUR VISIT

Access www.morelli.com.br
Information: 0800-0121455

Where to buy **our products**

Distributors

We have a network of more than 300 points of sale throughout the country.

Use the QR Code to find the one nearest you

Online Shop

On the Morelli website, we keep our portfolio available, with more than 1,700 items ready for delivery.

Our service channels

WhatsApp (15) **3238 8200**

0800 012 1455
0800 703 1455

sac@morelli.com.br

Opening hours from 8 a.m. to 5:48 p.m., Monday to Friday (except public holidays).
We have technical support trained to solve all your questions.

bluedental

Transform your clinic's daily routine. No wasting time.

A complete and dedicated manager for dentistry, BlueDental is the ideal partner for your practice.



Dr. Simone Cesar
@dentistamusical

The only complete, unlimited and 100% free of charge*

*In the available modules: Agenda, Patients, Finance and Stock


More than 130,000 registered patients in acting professional accounts

More than 350,000 clinics using resources for optimized management

Register and try

@bluedentalsistema

Access from anywhere, at any time and with any device connected to the Internet



Patient files

Manage and file patient records securely, creating personalized odontograms and storing photos and exams.

Scheduling

Make appointments and return visits, monitoring your schedule from anywhere with Internet access.

Financial management

Control financial and commissioning transactions digitally and generate reports easily.

Stock management

Manage your inputs and resources stock.

Sending messages

Via email or WhatsApp, it is possible to confirm, consult, congratulate the patient on their birthday and schedule a return visit.

Document automation

Receipts, contracts, odontograms, payment flows, medical records: all generated in just a few clicks.

Your data secured

In addition to maintaining a routine for backing up your data, it is encrypted using an SSL certificate. This ensures that no one else has access to the information registered in the system.

MORELLI CATALOG
Edition 12 • 2025

MORELLI



(15) 3238.8200 WHATSAPP
0800 0121455 • 0800 7031455
VENDAS@MORELLI.COM.BR
WWW.MORELLI.COM.BR

morelliortodontia
morelli.ortodontia
dentalmorelli
dental-morelli-ltda